



FILE 'CAPLUS' ENTERED AT 14:57:44 ON 13 AUG 2002  
THIS IS SUBJECT TO THE TERMS OF YOUR STN CUSTOMER AGREEMENT.  
PLEASE SEE "HELP USAGETERMS" FOR DETAILS.  
COPYRIGHT (C) 2002 AMERICAN CHEMICAL SOCIETY (ACS)

Copyright in the articles to which records in this database refer is held by the publishers listed in the PUBLISHER (PB) field (available for records published or updated in Chemical Abstracts after December 31, 1996), unless otherwise indicated in the original publications. The CA Lexicon is the copyrighted intellectual property of the American Chemical Society and is provided to assist you in searching databases on STN. Any dissemination, distribution, copying, or storing of this information, without the prior written consent of CAS, is strictly prohibited.

FILE COVERS: 1967 - 13 Aug 2002 VOL 137 ISS 7  
FILE LAST UPDATED: 12 Aug 2002 (20020812/ED)

This file contains CAS Registry Numbers for easy and accurate substance identification.

URL roles have been modified effective December 16, 2001. Please check your SDI profiles to see if they need to be revised. For information on CAS roles, enter HELP ROLES at an arrow prompt or use the CAS Roles thesaurus (URL field) in this file.

LN                    STF  
L10                  SCF 2026 AND 1006  
L11                  8066 SEA FILE=REGISTRY SSS FUL L8 AND L10  
L12                  11401 SEA FILE=CAPLUS ABB=ON L12  
L13                  2317 SEA FILE=CAPLUS ABB=ON SOLID SUPPORT# OBI  
L14                  6866 SEA FILE=CAPLUS ABB=ON MICROARRAY?/OBI OR MICRO(L)ARRAY?/OBI  
L15                  20 SEA FILE=CAPLUS ABB=ON L13(L) L15  
L16                  17 SEA FILE=CAPLUS ABB=ON L13(L) L17  
L17                  4 SEA FILE=CAPLUS ABB=ON (L19 AND L17) OR (L10 AND L15)

LN                    STF  
L10                  SCF 2026 AND 1006  
L11                  8066 SEA FILE=REGISTRY SSS FUL L8 AND L10  
L12                  11401 SEA FILE=CAPLUS ABB=ON L12  
L13                  2317 SEA FILE=CAPLUS ABB=ON SOLID SUPPORT# OBI  
L14                  6866 SEA FILE=CAPLUS ABB=ON MICROARRAY?/OBI OR MICRO(L)ARRAY?/OBI  
L15                  20 SEA FILE=CAPLUS ABB=ON L13(L) L15  
L16                  17 SEA FILE=CAPLUS ABB=ON L13(L) L17  
L17                  12 SEA FILE=CAPLUS ABB=ON (L19 OR L20) AND 9/SC, SX

*Section 9 =*

*Biotechnical Methods*

LN                    STF  
L10                  SCF 2026 AND 1006  
L11                  8066 SEA FILE=REGISTRY SSS FUL L8 AND L10  
L12                  11401 SEA FILE=CAPLUS ABB=ON L12  
L13                  2317 SEA FILE=CAPLUS ABB=ON SOLID SUPPORT# OBI  
L14                  6866 SEA FILE=CAPLUS ABB=ON MICROARRAY?/OBI OR MICRO(L)ARRAY?/OBI  
L15                  737 SEA FILE=CAPLUS ABB=ON L13(L)ANST/RL *Rel. AAST = analytical study*  
L16                  12 SEA FILE=CAPLUS ABB=ON (L15 OR L17) AND L20

L8 STF  
 L10 SCF 2026 AND 1006  
 L11 8068 SEA FILE=REGISTRY SSS FUL L8 AND L10  
 L12 21401 SEA FILE=CAPLUS ABB=ON L12  
 L13 737 SEA FILE=CAPLUS ABB=ON L12(L)ANST PL  
 L14 249468 SEA FILE=CAPLUS ABB=ON MICROF2/OBI  
 L15 135099 SEA FILE=CAPLUS ABB=ON ENA-OLD CT  
 L16 144371 SEA FILE=CAPLUS ABB=ON ENA-OLD CT  
 L17 6077 SEA FILE=CAPLUS ABB=ON IDENTIFIER CT  
 L18 6077 SEA FILE=CAPLUS ABB=ON POLYSACCHARIDES+OLD/CT  
 L19 131181 SEA FILE=CAPLUS ABB=ON L11111+OLD CT  
 L20 13416 SEA FILE=CAPLUS ABB=ON CVALENTS/OBI  
 L21 5580 SEA FILE=CAPLUS ABB=ON L11 OF L21 OR L33 OR L34 OR L35) (L) (L2  
 7 OR L36)  
 L22 9 SEA FILE=CAPLUS ABB=ON L22 AND L23

L8 STF  
 L10 SCF 2026 AND 1006  
 L11 8068 SEA FILE=REGISTRY SSS FUL L8 AND L10  
 L12 21401 SEA FILE=CAPLUS ABB=ON L12  
 L13 737 SEA FILE=CAPLUS ABB=ON MICRO SUPP EST#/OBI  
 L14 26610 SEA FILE=CAPLUS ABB=ON MICROFILIT OBI  
 L15 6066 SEA FILE=CAPLUS ABB=ON MICROARRAY/OBI OR MICRO(L)ARRAY?/OBI  
 L16 541 SEA FILE=CAPLUS ABB=ON L11(L) (L13 OR L16 OR L17)  
 L17 737 SEA FILE=CAPLUS ABB=ON L11(L)ANST PL  
 L18 1 SEA FILE=REGISTRY ABB=ON CYTOSINE/ON  
 L19 1 SEA FILE=REGISTRY ABB=ON GUANINE/ON  
 L20 10251 SEA FILE=CAPLUS ABB=ON L18 OF CYTOSINE/OBI  
 L21 37497 SEA FILE=CAPLUS ABB=ON L18 OF GUANINE/OBI  
 L22 3 SEA FILE=CAPLUS ABB=ON L18 OF L18) AND (L41 OR L42)

= s 123 or 124 or 126 or 138 or 143

L43 41 L23 OR L24 OF L11 OR L33 OR L41

= d ibib ans histstr 145 1-41; all item

L43 ANSWER 1 OF 41 CAPLUS COPYRIGHT 2001 A.M.

ABSESSION NUMBER: 1002:466034 CAPLUS

DOCUMENT NUMBER: 137:43915

TITLE: Method of attaching a biopolymer to a solid support  
using bromoacetamidoxilanes to functionalize the  
support

INVENTOR(S): Pirrung, Michael S.; Edenbaugh, Amy L.; Connors,  
Richard V.; Worden, John D.

PATENT ASSIGNEE(S): USA

SOURCE: U.S. Pat. Appl. Publ., 13 pp.

COLEN: USXKDD

DOCUMENT TYPE: Patent

LANGUAGE: English

FAMILY ACC. NUM. COUNT: 1

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
US 2002076832	A1	20010610	US 2001-371691	20010604
PRIORITY APPLN. INFO.:			US 2000-208493P	P 20000602
OTHER SOURCE(S):		MAFFAT 137:43915		

AB The present invention relates, in general, to a method of attaching a biopolymer to a solid support and, in particular, to a method of attaching a nucleic acid to a glass surface, and to reagents suitable for use in such a method. The invention further relates to the product produced by the present method and to kits comprising same. Clean microscope slides were silanized with N-(3-diethoxymethylsilylpropyl)bromoacetamide (prepn. given). Four oligonucleotides differing in only the nucleotide at their (free) 3'-ends were arrayed. When the array was treated with polymerase and fluoresceinated terminator, specific labeling of only the primer with perfect complementarity to the template was obsd.

13 3179-76-8, (3-Aminopropyl)methyldiethoxysilane 18306-79-1  
, 3-Aminopropyldimethylethoxysilane

RL: RCT (Reactant); RACT (Reactant or reagent)

(method of attaching biopolymers to **solid supports**

using bromoacetamidossilanes to functionalize supports)

EN 3179-76-8 CAPLUS

IN 1-Propanamine, 3-(diethoxymethylsilyl)- (9CI) (CA INDEX NAME)

OE:

Me Si (CH<sub>2</sub>)<sub>3</sub> NH<sub>2</sub>

OE:

EN 18306-79-1 CAPLUS

IN 1-Propanamine, 3-(ethoxydimethylsilyl)- (9CI) (CA INDEX NAME)

OE:

Me Si (CH<sub>2</sub>)<sub>3</sub> NH<sub>2</sub>

Me

EN 256352-86-0P 256352-87-1P 256352-89-3P  
437610-24-7P

RL: RCT (Reactant); SPN (Synthetic preparation); PREP (Preparation); RACT  
(Reactant or reagent)

(method of attaching biopolymers to **solid supports**

using bromoacetamidossilanes to functionalize supports)

EN 256352-86-0 CAPLUS

IN Acetamide, 2-bromo-N-[3-(diethoxymethylsilyl)propyl]- (9CI) (CA INDEX  
NAME)

OE: O

Me Si (CH<sub>2</sub>)<sub>3</sub> NH C CH<sub>2</sub>Br

OE:

EN 256352-87-1 CAPLUS

IN Acetamide, 2-bromo-N-[3-(ethoxydimethylsilyl)propyl]- (9CI) (CA INDEX  
NAME)

OEt O

Me Si (CH<sub>2</sub>)<sub>2</sub> NH C CH<sub>2</sub>Br

Me

EN 256350-44-3 CAPLUS

CN 1-Butanamine, 4-[methoxybis(1-methylethyl)silyl]- (9CI) (CA INDEX NAME)

OMe

i-Pr Si (CH<sub>2</sub>)<sub>4</sub> NH<sub>2</sub>

i-Pr

EN 447610-14-7 CAPLUS

CN Acetamide, 2-bromo-N-[4-(methoxybis(1-methylethyl)silyl]butyl]- (9CI) (CA INDEX NAME)

OMe O

i-Pr Si (CH<sub>2</sub>)<sub>4</sub> NH C CH<sub>2</sub>Br

i-Pr

L45 ANSWER 1 OF 41 CAPLUS COPYRIGHT 2002 ACS

ACCESSION NUMBER: 2002:173444 CAPLUS

DOCUMENT NUMBER: 136:129021

TITLE: High-density functional slide for biomolecule immobilization and preparation method thereof for use in high-efficiency bio-chip/microarray

INVENTOR(S): Ho, Chih-wei; Chow, Zu-sho; Fan, Bor-liuan; Tsao, Chia-huey; Pan, Chao-chi; Kuo, Wen-hsun; Chang, Yac-sun; Wu, Cheng-tao; Liu, Lu-ching

PATENT ASSIGNEE(S): Taiwan

SOURCE: U.S. Pat. Appl. Publ., 14 pp.

CODEN: SEXKNC

DOCUMENT TYPE: Patent

LANGUAGE: English

FAMILY ACC. NUM. COUNT: 1

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
US 2002018506	A1	2002/08/01	US 2001-886302	2001/04/18
PRIORITY APPLN. INFO.:			TW 2000-39118079 A 2000/09/04	

AB The invention features a method for prepg. a high-d. functional slide with ultra-thin layer by coating a sol-gel contg. amine-group bearing silanes and a soln. contg. polyaldehyde groups onto an org. or inorg. substrate, resp. The resulting slide is useful in the prepn. of highly homogeneous functional-group slides and the high-d. and high-efficiency bio-chip/microarray. In one preferred embodiment of the present invention, the polyaldehyde polymer is prepd. via the graft co-polymn. of polyvinylalc.-based polyaldehyde. Therefore, the present invention also provides a polyvinylalc.-based polyaldehyde graft copolymer, which is prepd. by the following steps: (a) dissolving polyvinylalc. in water to

form a polymeric soln.; (b) adding the monomer of allyl alc. and acrolein to the polymeric soln. under anaerobic conditions; and (c) adding ceric ammonium nitrate to the soln. for catalysis. The polyvinylalc.-based polyaldehyde graft copolymer comprises 2-10 (w/w) polyvinylalc., 2-10 (vol. vol.) monomer of acrolein and 1-5 (vol./vol.) monomer of allyl alc.

IT 919-30-2, Aminopropyltriethoxysilane

RL: DEV (Device component use); USES (Uses)

(APTES, sol-gel; high-d. functional slide for biomol. immobilization and prepn. method thereof for high-efficiency bio-chips; **microarray**)

RN 919-30-2 CAPLUS

CN 1-Propanamine, 3-(triethoxysilyl)- (9CI) (CA INDEX NAME)

OEt

EO Si (CH<sub>2</sub>)<sub>3</sub> NH<sub>2</sub>

OEt

L4: ANSWER 3 OF 41 CAPLUS COPYRIGHT 2002 ACS

ACCESSION NUMBER: 2002:90792 CAPLUS

DOCUMENT NUMBER: 136:275611

TITLE: Characteristics of DNA **microarrays**

fabricated on various aminosilane layers

AUTHOR(S): Oh, Seon Jin; Cho, Sung Ju; Kim, Chang Ok; Park, Joon Won

CORPORATE SOURCE: Center for Integrated Molecular Systems, Department of Chemistry, Division of Molecular and Life Sciences, Pohang University of Science and Technology, Pohang, 790-784, S. Korea

SOURCE: Langmuir (2002), 18(1), 1764-1769

CODEN: LANGDE; ISSN: 0743-7463

PUBLISHER: American Chemical Society

DOCUMENT TYPE: Journal

LANGUAGE: English

AB Four kinds of aminosilane layers on glass slides or silicon wafers were prepd. The amine densities of the layers prepd. with aminopropyltriethoxymethylsilane (APDES), aminopropylmonoethoxydimethylsilane (APMES), a mixt. of (aminopropyl)triethoxysilane (APTES) and n-butyltrimethoxysilane (n-BTMO) (vol./vol. = 1:10) are 4.0(±0.8), 1.3(±0.1), and 0.30(±0.6) amine/nm<sup>2</sup>, resp. A substrate with much higher amine d., i.e., 40(±1) amines/nm<sup>2</sup> was also prepd. by allowing aziridine to polymerize on the APDES-treated substrate. AFM revealed that APDES-, APMES-, and APTES/n-BTMO-treated surfaces were relatively flat; on the other hand, an aziridine-treated surface showed embossed morphol. The amine substrates were allowed to react with a heterobifunctional linker succinimidyl 4-maleimido butyrate (SMB), and subsequently pentadecadecoxynucleotides were microarrayed on the SMB-treated substrates. Characteristics of the DNA microarrays including the dynamic range, the mismatch discrimination efficiency, and so forth were examd. Noteworthy, DNA microarrays on the aziridine-polymd. substrate showed much higher fluorescence intensity. At the same time, DNA microarrays from these four substrates were able to discriminate internal- and terminal-mismatched pairs, but the fluorescence ratio was far from the one that thermodyn. implies.

IT 919-30-2, APTES 3179-76-8 18306-79-1

RL: ANS (Analytical role, unclassified); PEP (Physical, engineering or chemical process); PTF (Physical process); ANST (Analytical study); PROC (Process)

(DNA microarrays fabricated on various aminosilane layers)

RN 319-30-2 CAPLUS  
CN 1-Propanamine, 3-(triethoxysilyl)- (9CI) (CA INDEX NAME)

OEt

EtO Si (CH<sub>3</sub>)<sub>3</sub> NH<sub>2</sub>

OEt

RN 3179-76-8 CAPLUS  
CN 1-Propanamine, 3-(diethoxymethylsilyl)- (9CI) (CA INDEX NAME)

OEt

Me Si (CH<sub>3</sub>)<sub>3</sub> NH<sub>2</sub>

OEt

RN 18390-79-1 CAPLUS  
CN 1-Propanamine, 3-(ethoxydimethylsilyl)- (9CI) (CA INDEX NAME)

OEt

Me Si (CH<sub>3</sub>)<sub>3</sub> NH<sub>2</sub>

Me

REFERENCE COUNT: 45 THERE ARE 45 CITED REFERENCES AVAILABLE FOR THIS  
RECORD. ALL CITATIONS AVAILABLE IN THE RE FORMAT

145 ANSWER 4 OF 41 CAPLUS COPYRIGHT 2002 ACS

ACCESSION NUMBER: 2002:51931 CAPLUS

DOCUMENT NUMBER: 136:80850

TITLE: Compositions and methods for array-based genomic  
nucleic acid analysis of biological molecules

INVENTOR(S): Bradley, Allan; Cai, Wei-Wen; Marathi, Upendra

PATENT ASSIGNEE(S): UK

SOURCE: U.S. Pat. Appl. Publ., 25 pp., Cont.-in-part of U.S.  
Ser. No. 546,035.

CODEN: USXKCO

DOCUMENT TYPE: Patent

LANGUAGE: English

FAMILY APP. NUM. COUNT: 2

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
US 2002006623	A1	20020117	US 2001-853343	20010510
US 6048695	A	20000411	US 1998-7187	19980504
PRIORITY APPLN. INFO.:			US 1999-71876	A2 19980504
			US 2000-540085	A2 20000410

OTHER SOURCE(S): WAFPAT 136:80850

AB The invention provides biomol. modified by reaction with a compd.  
having the formula: R1-X-R2, wherein R1 is a cyclic ether group or an  
amino group, R2 is an alkoxy-silane group and X is a moiety chem. suitable  
for linking the cyclic ether group or the amino group to the alkoxy-silane

group. The invention also provides arrays, or "biochips," comprising these modified biol. mols. Also provided are methods for making and using these compns.

919-30-2, 3-Aminopropyltriethoxysilane 2530-83-8,

3-Glycidyloxypropyltrimethoxysilane

RL: ANS (Analytical reagent use); BIO (Biological use, unclassified);

ANST (Analytical study); BIOL (Biological study); USES (Uses)

compr. and methods for array-based genomic nucleic acid anal. or biol. mols.

919-30-2 CAPLUS

3-Propylamine, 3-(triethoxysilyl)- (9CI) (CA INDEX NAME)

DEt

919-30-2 (CH<sub>2</sub>)<sub>3</sub> NH<sub>2</sub>

DEt

919-30-2 CAPLUS

3-Silane, trimethoxy[3-(oxiranylmethoxy)propyl]- (9CI) (CA INDEX NAME)

OMe

CH<sub>2</sub> O (CH<sub>2</sub>)<sub>3</sub> Si OMe

OMe

L45 ANSWER 5 OF 41 CAPLUS COPYRIGHT 2012 ACS

ACCESSION NUMBER: 2002:31489 CAPLUS

DOCUMENT NUMBER: 156:31739

TITLE: Improved combination of microporous membrane and solid support for micro-analytical diagnostic applications

PATENT ASSIGNEE(S): Cuno, Inc., USA

SOURCE: ECT Int. Appl., 39 pp.

CODEN: PIXXD2

DOCUMENT TYPE: Patent

LANGUAGE: English

FAMILY ACC. NUM. COUNT: 2

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
WO 2002004477	A2	20020117	WO 2001-0321210	20010603
WO 2002004477	C1	20020629		

W: AU, BR, JP

RW: AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE, TR

US 2002086307 A1 20020704 US 2001-393102 20010703

PRIORITY APPLN. INFO.: US 2000-216390P P 20000706

AB The invention concerns an improved combination microporous membrane and solid support for use in micro-anal. diagnostic applications is disclosed. Specifically, a multi-cell substrate useful for carrying a microarray of biol. polymers on the surface thereof including a multi-cell substrate having a porous membrane formed by a phase inversion process effectively attached by covalent bonding through a surface treatment to a substrate that preps. the substrate to sufficiently, covalently bond to the microporous membrane formed by a phase inversion process such that the combination produced thereby is useful in microarray applications and



wherein the porous nylon multi-cell substrate is covalently bonded to a solid base member, such as, for example, a glass or Mylar microscope slide, such that the combination produced thereby is useful in microarray applications. App. for fabricating a multi-cell substrate is also disclosed. Diagrams describing the app. are given.

IT 919-30-2, 3-Aminopropyltriethoxysilane 1760-24-3,  
N-(2-Aminoethyl)-3-aminopropyltrimethoxysilane 2530-83-8,  
3-Glycidyloxypropyltrimethoxysilane  
PL: NUU (Other use, unclassified); USES (Uses)  
(Improved combination of microporous membrane and **solid support** for micro-anal. diagnostic applications)  
RN 919-30-2 CAPLUS  
CN 1-Propanamine, 3-(triethoxysilyl)- (9CI) (CA INDEX NAME)

OMe

EtO Si (CH<sub>2</sub>)<sub>3</sub> NH<sub>2</sub>

OMe

RN 1760-24-3 CAPLUS  
CN 1,2-Ethanediamine, N-[3-(trimethoxysilyl)propyl]- (9CI) (CA INDEX NAME)

OMe

MeO Si (CH<sub>2</sub>)<sub>3</sub> NH CH<sub>2</sub> CH<sub>2</sub> NH<sub>2</sub>

OMe

RN 2530-83-8 CAPLUS  
CN Silane, trimethoxy[3-(oxiranylmethoxy)propyl]- (9CI) (CA INDEX NAME)

OMe

CH<sub>2</sub> O (CH<sub>2</sub>)<sub>3</sub> Si OMe

OMe

L45 ANSWER 6 OF 41 CAPLUS COPYRIGHT 2002 ACS  
ACCESSION NUMBER: 1801:090855 CAPLUS  
DOCUMENT NUMBER: 135:161917  
TITLE: An Efficient Binding Chemistry for Glass Polynucleotide Microarrays  
AUTHOR(S): Lee, Paul H.; Sawa, Samuel P.; Mcdrusan, Zora; Arnold, Lyle T., Jr.; Reynolds, Mark A.  
CORRESPONDENCE SOURCE: Incyte Genomics, Microarray Research and Development, Fremont, CA, 94555, USA  
SOURCE: Bioconjugate Chemistry (2002), 13(1), 97-103  
CODEN: BCCHEH; ISSN: 1043-1802  
PUBLISHER: American Chemical Society  
DOCUMENT TYPE: Journal  
LANGUAGE: English

AB A variety of methods have been described for making synthetic polynucleotide microarrays. These include in situ synthesis directly on the array surface, for example, by photolithog. or ink-jet printing technologies, and the application of presynthesized polynucleotides that

are derivatized with various nucleophiles or electrophiles. In the latter case, a variety of surface chemistries have been developed, and several are available com. These chemistries must be compatible with nanoliter-scale vols. of polynucleotide reagents, which contact the array over a small portion of their surface. We reasoned that a three-dimensional polymer coating could potentially offer greater surface contact and higher binding efficiency. Here we describe a polyethylenimine-based coating chem. that provides exceptional binding and hybridization characteristics. In our preferred process, size-fractionated polyethylenimine polymers are cross-linked onto an aminopropylsilanated glass surface in the presence of cyanuric chloride. The resulting three-dimensional coating binds polynucleotides through a mixt. of covalent and noncovalent interactions as evidenced by comparisons between 5'-aminoalkyl modified and unmodified polynucleotides. Binding and hybridization comparisons are presented including analogous two-dimensional electrophilic and electrostatic chemistries.

BT 13822-56-5, 3-Aminopropyltrimethoxysilane

RE: RCT (Reactant); RACT (Reactant or reagent)

efficient binding chem. for glass polynucleotide **microarrays**, synthesis and characterization of glass surface coatings)

BT 13822-56-5 CAPLUS

BT 1-Propanamine, 3-(trimethoxysilyl)- (PCI) (CA INDEX NAME)

OMe

MeO Si (CH<sub>2</sub>)<sub>3</sub> NH<sub>2</sub>

OMe

REFERENCE COUNT: 17 THERE ARE 17 CITED REFERENCES AVAILABLE FOR THIS RECORD. ALL CITATIONS AVAILABLE IN THE RE FORMAT

BT ANSWER 7 OF 41 CAPLUS COPYRIGHT 2002 ACS

ACCESSION NUMBER: 2001:363771 CAPLUS

DOCUMENT NUMBER: 136:163471

TITLE: HPLC of some nucleosides and bases on p-tert-butyl-calix[6]arene-bonded silica gel stationary phase

AUTHOR(S): Xiao, Yu-Kun; Xiao, Xiang-Zhu; Feng, Yu-Qi; Wang, Zhong-Hua; Da, Shi-Lu

CORPORATE SOURCE: College of Life Sciences and Department of Chemistry, Wuhan University, Wuhan, 430072, Peop. Rep. China

SOURCE: Journal of Liquid Chromatography & Related Technologies (2001), 24(19), 2923-2942

CODEN: JLCSTC; ISSN: 1082-6076

PUBLISHER: Marcel Dekker, Inc.

DOCUMENT TYPE: Journal

LANGUAGE: English

AB The high-performance liq. chromatog. behavior of some nucleosides and bases was studied on a new p-tert-butyl-calix[6]arene-bonded silica gel stationary phase. The effect of mobile phase variables, such as ionic strength, methanol content, and pH on their chromatog. behavior was investigated. Some nucleosides and bases were successfully sepd. on the new stationary phase. Their retention behavior was compared with that on both Porbax C18 phase and-(6-ethylenediamine)propyl-triethoxysilane-bonded silica gel. The results indicate that the new stationary phase behaves as a reversed-phase packing, but its hydrophobicity is much weaker than that of Porbax C18 phase. The retention mechanism on the new stationary phase was also discussed.

BT 5089-72-5D, reaction products with polycalixarene acetic acid chlorides

PL: AFU (Analytical role, unclassified); ANST (Analytical study)  
(HPLC of nucleosides and bases on p-tert-butyl-calix[6]arene-bonded  
silica gel stationary phase)

EN 5989-72-5 CAPLUS  
CN 1,2-Ethanediamine, N-[3-(triethoxysilyl)propyl]- (9CI) (CA INDEX NAME)

OE:

EtO Si (CH<sub>3</sub>)<sub>3</sub> NH CH<sub>2</sub> CH<sub>2</sub> NH<sub>2</sub>

OE:

IT 71-30-7, Cytosine 73-40-5, Guanine

PL: PEI (Physical, engineering or chemical process); PYP (Physical  
process;; FPOC (Process)  
(HPLC of nucleosides and bases on p-tert-butyl-calix[6]arene-bonded  
silica gel stationary phase)

EN 71-30-7 CAPLUS  
CN (1H)-Pyrimidinone, 4-amino- (9CI) (CA INDEX NAME)

C     H  
      N     NH<sub>2</sub>

N

EN 73-40-5 CAPLUS  
CN 6H-Purin-6-one, 2-amino-1,7-dihydro- (9CI) (CA INDEX NAME)

H<sub>2</sub>N     H     N  
          N     NH

O

REFERENCE COUNT: 26     THERE ARE 26 CITED REFERENCES AVAILABLE FOR THIS  
RECORD. ALL CITATIONS AVAILABLE IN THE RE FORMAT

L45 ANSWER 8 OF 41 CAPLUS CCFYF1GHT 2002 ACS  
ACCESSION NUMBER: 2001:748054 CAPLUS  
DOCUMENT NUMBER: 135:299485  
TITLE: Compositions and methods for detecting and quantifying  
gene expression in microarrays  
INVENTOR(S): Lowe, David G.; Marsters, James C., Jr.; Robbie,  
Edward P.; Smith, Victoria  
PATENT ASSIGNEE(S): Genentech, Inc., USA  
SOURCE: PCT Int. Appl., 54 pp.  
CODEN: PIMXD2  
DOCUMENT TYPE: Patent  
LANGUAGE: English  
FAMILY ACC. NUM. COUNT: 1  
PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
-----	----	-----	-----	-----

WI: 1076160 A2: 20011911 WO: 2001-031492 2001/001  
WI: 1076160 A3: 20023501

WI: AE, AG, AL, AM, AT, AU, AZ, BA, BB, BC, BD, BE, BF, BG, BH, BI, BJ, BK, BL, BM, BN, BO, BP, BQ, BR, BS, BT, BU, BV, BW, BX, BY, BZ, CA, CH, CN, CO, CP, CR, CU, CV, CW, CX, CY, CZ, DA, DB, DC, DD, DE, DF, DG, DH, DI, DJ, DK, DL, DM, DN, DO, DP, DQ, DR, DS, DT, DU, DV, DW, DX, DY, DZ, EA, EB, EC, ED, EE, EF, EG, EH, EI, EJ, EK, EL, EM, EN, EO, EP, EQ, ER, ES, ET, EU, EV, EW, EX, EY, EZ, FA, FB, FC, FD, FE, FF, FG, FH, FI, FJ, FK, FL, FM, FN, FO, FP, FQ, FR, FS, FT, FU, FV, FW, FX, FY, FZ, GA, GB, GC, GD, GE, GF, GH, GI, GJ, GK, GL, GM, GN, GO, GP, GQ, GR, GS, GT, GU, GV, GW, GX, GY, GZ, HA, HB, HC, HD, HE, HF, HG, HH, HI, HJ, HK, HL, HM, HN, HO, HP, HQ, HR, HS, HT, HU, HV, HW, HX, HY, HZ, IA, IB, IC, ID, IE, IF, IG, IH, II, IJ, IK, IL, IM, IN, IO, IP, IQ, IR, IS, IT, IU, IV, IW, IX, IY, IZ, JA, JB, JC, JD, JE, JF, JG, JH, JI, JJ, JK, JL, JM, JN, JO, JP, JQ, JR, JS, JT, JU, JV, JW, JX, JY, JZ, KA, KB, KC, KD, KE, KF, KG, KH, KI, KJ, KK, KL, KM, KN, KO, KP, KQ, KR, KS, KT, KU, KV, KW, KX, KY, KZ, LA, LB, LC, LD, LE, LF, LG, LH, LI, LJ, LK, LM, LN, LO, LP, LQ, LR, LS, LT, LU, LV, LW, LX, LY, LZ, MA, MB, MC, MD, ME, MF, MG, MH, MI, MJ, MK, ML, MM, MN, MO, MP, MQ, MR, MS, MT, MU, MV, MW, MX, MY, MZ, NA, NB, NC, ND, NE, NF, NG, NH, NI, NJ, NK, NL, NM, NO, NP, NQ, NR, NS, NT, NU, NV, NW, NX, NY, NZ, OA, OB, OC, OD, OE, OF, OG, OH, OI, OJ, OK, OL, OM, ON, OO, OP, OQ, OR, OS, OT, OU, OV, OW, OX, OY, OZ, PA, PB, PC, PD, PE, PF, PG, PH, PI, PJ, PK, PL, PM, PN, PO, PP, PQ, PR, PS, PT, PU, PV, PW, PX, PY, PZ, QA, QB, QC, QD, QE, QF, QG, QH, QI, QJ, QK, QL, QM, QN, QO, QP, QQ, QR, QS, QT, QU, QV, QW, QX, QY, QZ, RA, RB, RC, RD, RE, RF, RG, RH, RI, RJ, RK, RL, RM, RN, RO, RP, RQ, RR, RS, RT, RU, RV, RW, RX, RY, RZ, SA, SB, SC, SD, SE, SF, SG, SH, SI, SJ, SK, SL, SM, SN, SO, SP, SQ, SR, SS, ST, SU, SV, SW, SX, SY, SZ, TA, TB, TC, TD, TE, TF, TG, TH, TI, TJ, TK, TL, TM, TN, TO, TP, TQ, TR, TS, TU, TV, TW, TX, TY, TZ, UA, UB, UC, UD, UE, UF, UG, UH, UI, UJ, UK, UL, UM, UN, UO, UP, UQ, UR, US, UT, UV, UW, UX, UY, UZ, VA, VB, VC, VD, VE, VF, VG, VH, VI, VJ, VK, VL, VM, VN, VO, VP, VQ, VR, VS, VT, VU, VV, VW, VX, VY, VZ, WA, WB, WC, WD, WE, WF, WG, WH, WI, WJ, WK, WL, WM, WN, WO, WP, WQ, WR, WS, WT, WU, WV, WW, WX, WY, WZ, XA, XB, XC, XD, XE, XF, XG, XH, XI, XJ, XK, XL, XM, XN, XO, XP, XQ, XR, XS, XT, XU, XV, XW, XX, XY, XZ, YA, YB, YC, YD, YE, YF, YG, YH, YI, YJ, YK, YL, YM, YN, YO, YP, YQ, YR, YS, YT, YU, YV, YW, YX, YY, YZ, ZA, ZB, ZC, ZD, ZE, ZF, ZG, ZH, ZI, ZJ, ZK, ZL, ZM, ZN, ZO, ZP, ZQ, ZR, ZS, ZT, ZU, ZV, ZW, ZX, ZY, ZZ.

WI: GH, GI, GJ, GK, GL, GM, GN, GO, GP, GQ, GR, GS, GT, GU, GV, GW, GX, GY, GZ, HA, HB, HC, HD, HE, HF, HG, HH, HI, HJ, HK, HL, HM, HN, HO, HP, HQ, HR, HS, HT, HU, HV, HW, HX, HY, HZ, IA, IB, IC, ID, IE, IF, IG, IH, II, IJ, IK, IL, IM, IN, IO, IP, IQ, IR, IS, IT, IU, IV, IW, IX, IY, IZ, JA, JB, JC, JD, JE, JF, JG, JH, JI, JJ, JK, JL, JM, JN, JO, JP, JQ, JR, JS, JT, JU, JV, JW, JX, JY, JZ, KA, KB, KC, KD, KE, KF, KG, KH, KI, KJ, KK, KL, KM, KN, KO, KP, KQ, KR, KS, KT, KU, KV, KW, KX, KY, KZ, LA, LB, LC, LD, LE, LF, LG, LH, LI, LJ, LK, LM, LN, LO, LP, LQ, LR, LS, LT, LU, LV, LW, LX, LY, LZ, MA, MB, MC, MD, ME, MF, MG, MH, MI, MJ, MK, ML, MM, MN, MO, MP, MQ, MR, MS, MT, MU, MV, MW, MX, MY, MZ, NA, NB, NC, ND, NE, NF, NG, NH, NI, NJ, NK, NL, NM, NO, NP, NQ, NR, NS, NT, NU, NV, NW, NX, NY, NZ, OA, OB, OC, OD, OE, OF, OG, OH, OI, OJ, OK, OL, OM, ON, OO, OP, OQ, OR, OS, OT, OU, OV, OW, OX, OY, OZ, PA, PB, PC, PD, PE, PF, PG, PH, PI, PJ, PK, PL, PM, PN, PO, PP, PQ, PR, PS, PT, PU, PV, PW, PX, PY, PZ, QA, QB, QC, QD, QE, QF, QG, QH, QI, QJ, QK, QL, QM, QN, QO, QP, QQ, QR, QS, QT, QU, QV, QW, QX, QY, QZ, RA, RB, RC, RD, RE, RF, RG, RH, RI, RJ, RK, RL, RM, RN, RO, RP, RQ, RR, RS, RT, RU, RV, RW, RX, RY, RZ, SA, SB, SC, SD, SE, SF, SG, SH, SI, SJ, SK, SL, SM, SN, SO, SP, SQ, SR, SS, ST, SU, SV, SW, SX, SY, SZ, TA, TB, TC, TD, TE, TF, TG, TH, TI, TJ, TK, TL, TM, TN, TO, TP, TQ, TR, TS, TU, TV, TW, TX, TY, TZ, UA, UB, UC, UD, UE, UF, UG, UH, UI, UJ, UK, UL, UM, UN, UO, UP, UQ, UR, US, UT, UV, UW, UX, UY, UZ, VA, VB, VC, VD, VE, VF, VG, VH, VI, VJ, VK, VL, VM, VN, VO, VP, VQ, VR, VS, VT, VU, VV, VW, VX, VY, VZ, WA, WB, WC, WD, WE, WF, WG, WH, WI, WJ, WK, WL, WM, WN, WO, WP, WQ, WR, WS, WT, WU, WV, WW, WX, WY, WZ, XA, XB, XC, XD, XE, XF, XG, XH, XI, XJ, XK, XL, XM, XN, XO, XP, XQ, XR, XS, XT, XU, XV, XW, XX, XY, XZ, YA, YB, YC, YD, YE, YF, YG, YH, YI, YJ, YK, YL, YM, YN, YO, YP, YQ, YR, YS, YT, YU, YV, YW, YX, YY, YZ, ZA, ZB, ZC, ZD, ZE, ZF, ZG, ZH, ZI, ZJ, ZK, ZL, ZM, ZN, ZO, ZP, ZQ, ZR, ZS, ZT, ZU, ZV, ZW, ZX, ZY, ZZ.

US 2002081597 A1 20020627 US 2001-823649 20010830

PRIORITY APPLN. INFO.:

US 2000-1937-07P P 20000831

AB: Comps. and methods for improving detection sensitivity in nucleic acid microarray anal. are disclosed, including methods of purifying nucleic acids, methods of synthesizing fluorescent DNA probes, methods of hybridization, and methods of activating a substrate for target mol. attachment. The comps. and methods of this invention include synthesis of cDNA, sDNA, or cRNA probes from cellular RNA by in vitro transcription and/or a single-round of reverse transcription with incorporation of fluorochromes. Specific procedures for microarray slide prepn. to decrease background fluorescence are given. For example, silanization of glass slides with toluene as the solvent is preferred. In addn., unmodified polynucleotides can attach to a glass slide treated with 3-aminopropyltriethoxysilane followed by phenylene diisothiocyanate. Modified target DNA can also be synthesized using PCR primers which contain a primary amine and an alkyl linker attached to the 5'-end. The modified target DNA is then reacted with activated silanized glass slides. Microarray hybridization buffers contg. alkylammonium salts, dimethylsulfoxide and formamide and lacking the detergent sodium dodecyl sulfate also improved the detection sensitivity. The invention is illustrated with microarrays hybridized with fluorescent probes synthesized from very small quantities of RNA isolated from microdissected tumor cells, paraffin-embedded liver and colon tissue, fresh frozen liver tissue, and fresh frozen colon tissue. The microarray expts. were designed to compare tissue sample prepn. methods and gene expression in tumor vs. healthy tissues. An example of the sensitivity of these methods shows a microarray hybridized with sDNA probes from one round of amplification of 2 pg of RNA from an ovarian carcinoma cell line.

IT: 919-30-2, 3-Aminopropyltriethoxysilane

EL: BUN (Biological use, unclassified); DEV (Device component use); RCT (Reactant); BIOL (Biological study); RACT (Reactant or reagent); USES (Uses)

(comps. and methods for detecting and quantifying gene expression in microarrays)

EN: 919-30-2 CAPLUS

IN: 1-Propanamine, 3-(triethoxysilyl)- (9CI) (CA INDEX NAME)

OE:

SI: (CH<sub>3</sub>)<sub>3</sub> NH<sub>2</sub>

OF:

LI: ANSWER 9 OF 41 CAPLUS COPYRIGHT 2001 ACS

ACCESSION NUMBER: 2001:718245 CAPLUS

COMMENT NUMBER: 133:26:66

TITLE: Polymer coated surfaces for microarray applications

INVENTOR(S): Arnold, Lyne J., Jr.; Sawan, Samuel P.; Lee, Paul H.

PATENT ASSIGNEE(S): Incyte Pharmaceuticals, Inc., USA

SOURCE: PCT Int. Appl., 2# pp.  
 CODEN: PIXXD2  
 DOCUMENT TYPE: Patent  
 LANGUAGE: English  
 FAMILY ACC. NUM. COUNT: 1  
 PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
WO 1001070641	A1	10010627	WO 2001-038933	10010321
W:	AE, AG, AL, AM, AT, AU, AZ, BA, BB, BG, BP, BY, BS, CA, CH, CN, CO, CP, CU, CC, DE, DK, EM, EZ, EE, ES, FI, GB, GR, GE, GH, GM, HE, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, LC, LK, LR, LS, LT, LU, LV, MA, MI, MS, MP, MN, MW, MX, NL, NO, NI, PL, PT, RO, RU, SD, SE, SG, SI, SK, SL, TJ, TM, TR, TT, TZ, UA, UG, UZ, VN, TH, VA, ZW, AM, AN, AP, AR, AT, AU, BA, BB, BG, BR, BS, BY, BW, CA, CH, CY, CZ, DE, DK, ES, FI, FR, GB, GR, HE, IE, IT, LU, MA, NL, PT, SE, TR, BF, BJ, CF, CG, CI, CM, GA, GN, GW, ML, MR, NE, NI, PD, TG			
US 6418711	B1	100010702	US 2000-532413	20000322
US 100107500	A1	100010702	US 2001-775413	20010201
US 6387681	B2	100010514		

PRIORITY APPLN. INFO.: US 2000-532413 A 20000322  
 OTHER SOURCE(S): MARIAT 135:260660

AB Methods are provide for modifying a solid support, such as a glass slide, by silylating with an agent having the formula  $H_2N(CH_2)_nSiX_3$  ( $n = 1-10$ ,  $X =$  independently chosen from Me, Et, Pr, Bu, I), then activating with a crosslinking reagent, followed by reacting with an amine-contg. polymer. The support can optionally be reacted with a crosslinking reagent again. The support thus modified may be used to make arrays and microarrays where a plurality of targets are stably assod. with the support and arranged in a defined manner. Thus, glass slides were silylated with 3-aminopropyltrimethoxysilane. The silylated slides were reacted with cyanuric chloride then with PEI, polylysine, or polyhistidine. 5'-Aminoalkyl-terminated oligonucleotides were spotted on such slides and used in hybridization assays.

IT 13822-56-5, 3-Aminopropyltrimethoxysilane  
 EL: PCT (Reactant); PACT (Reactant or reagent)  
 (polymer coated surfaces for **microarray** applications)  
 FI: 13822-56-5 CAPLUS  
 CH 1-Propanamine, 3-(trimethoxysilyl)- (931) (CA INDEX NAME)

CM-

MeO Si (CH<sub>2</sub>)<sub>3</sub> NH<sub>2</sub>

CM-

REFERENCE COUNT: 1 THERE ARE 5 PCTE REFERENCES AVAILABLE FOR THIS PCTED. AND CITATIONS AVAILABLE IN THE RE FORMAT

L45 ANSWER 10 OF 41 CAPLUS COPYRIGHT 2001 ACS

ACCESSION NUMBER: 1001:038933 CAPLUS

DOCUMENT NUMBER: 10010321

TITLE: Dendrimer-activated **solid supports**  
 for nucleic acid and protein **microarrays**

AUTHOR(S): Benters, R.; Niemeyer, C. M.; Wehrle, D.

CORPORATE SOURCE: Institut of Organic and Macromolecular Chemistry,  
 University Bremen, Bremen, 28359, Germany

SOURCE: ChemBioChem (2001), 2(9), 686-694

CODEN: CBBHEX; ISSN: 1439-4227

Searched by Barb O'Bryen, STIC 1001-4291

PUBLISHER: Wiley-VCH Verlag GmbH  
 DOCUMENT TYPE: Journal  
 LANGUAGE: English

AB The generation of chem. activated glass surfaces is of increasing interest for the prodn. of microarrays contg. DNA, proteins, and low-mol.-wt. components. We here report on a novel surface chem. for highly efficient activation of glass slides. Our method is based on the initial modification of glass with primary amino groups using a protocol, specifically optimized for high aminosilylation yields, and in particular, for homogeneous surface coverages. In a following step the surface amino groups are activated with a nonbifunctional linker, such as disuccinimidylglutarate (DSG) or 1,4-phenylenedithiocyanate (PDITC), and then allowed to react with a starburst dendrimer that contains 64 primary amino groups in its outer sphere. Subsequently, the dendritic monomers are activated and crosslinked with a nonbifunctional spacer, either DSG or PDITC. This leads to the formation of a thin, chem. reactive polymer film, covalently affixed to the glass substrate, which can directly be used for the covalent attachment of amino-modified components, such as oligonucleotides. The resulting DNA microarrays were studied by means of nucleic acid hybridization expts. using fluorophore-labeled complementary oligonucleotide targets. The results indicate that the novel dendrimer-activated surfaces display a surface coverage with capture oligomers about twofold greater than that with conventional microarrays contg. linear chem. linkers. In addn., the expts. suggest that the hybridization occurs with decreased steric hindrance, likely a consequence of the long, flexible linker chain between the surface and the DNA oligomer. The surfaces were found to be resistant against repeated alk. regeneration procedures, which is likely a consequence of the crosslinked polymeric structure of the dendrimer film. The high stability allows multiple hybridization expts. without significant loss of signal intensity. The versatility of the dendrimer surfaces is also demonstrated by the covalent immobilization of streptavidin as a model protein.

IT 392661-75-5 392661-76-6


EL: ARU (Analytical role, unclassified); DEV (Device component use);

ANST (Analytical study); USES (Uses)

condensation on silica; dendrimer-activated solid supports for nucleic acid and protein microarrays)

BN 392661-75-5 CAPLUS

IN Pentanamide, 5-[(2,5-dioxo-1-pyrrolidinyloxy]-5-oxo-N-[3-(triethoxysilyl)propyl]- (9CI) (CA INDEX NAME)



BN 392661-76-6 CAPLUS

IN Thiourea, N-(4-isothiocyanatophenyl)-N'-[3-(triethoxysilyl)propyl]- (9CI) (CA INDEX NAME)

S OEt

NH C NH (CH<sub>2</sub>)<sub>3</sub> Si OEt

Et

S C N

REFERENCE COUNT: 34 THERE ARE 34 CITED REFERENCES AVAILABLE FOR THIS  
 RECORD. ALL CITATIONS AVAILABLE IN THE RE FORMAT

145 ANSWER 11 OF 41 CAPLUS COPYRIGHT 2002 ACS  
 ACCESSION NUMBER: 2001:611629 CAPLUS  
 DOCUMENT NUMBER: 133:17767  
 TITLE: Linear microarrays  
 INVENTOR(S): Helann, Timothy W.; Park, Sang Chul  
 PATENT ASSIGNEE(S): Incyte Genomics, Inc., USA  
 SOURCE: 3.77, 11 pp.  
 KEYW: USXNAM  
 DOCUMENT TYPE: Patent  
 LANGUAGE: English  
 FAMILY ACC. NUM. COUNT: 1  
 PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
US 6277613	B1	20010311	US 1998-165465	19981002
US 2001072065	A1	20020611	US 2001-935570	20010810

PRIORITY AFFILI. INFO.: US 1998-165465 A1 19931002

AB The present invention provides a method and a compn. for detecting the levels of a plurality of biomed. probes in a sample. In particular, the invention relates to a hybridization compn. for detecting the presence or levels of different polynucleotide sequences in a sample. A YP3 59mer labeled at the 3'-end with a Cy3 fluorescent dye was immobilized on epoxide-coated glass beads. A capillary tube was packed with the beads sepd. by alternating unmodified beads to prep. a glass bead array.

IT 2530-83-8, 3-Glycidioxypropyl-trimethoxysilane  
 FL: FCT (Reactant); FACT (Reagent or reagent)  
 (linear microarrays

FI 2530-83-3 CAPLUS

CU Silane, trimethoxy[3-(trimethoxypropyl)- (OCI) (CA INDEX NAME)

C OMe

CH<sub>2</sub> O (CH<sub>2</sub>)<sub>3</sub> Si OMe

OMe

REFERENCE COUNT: 4 THERE ARE 4 CITED REFERENCES AVAILABLE FOR THIS  
 RECORD. ALL CITATIONS AVAILABLE IN THE RE FORMAT

145 ANSWER 12 OF 41 CAPLUS COPYRIGHT 2002 ACS  
 ACCESSION NUMBER: 2001:57613 CAPLUS  
 DOCUMENT NUMBER: 133:194833  
 TITLE: Oligonucleotides form a duplex with non-helical properties on a positively charged surface  
 AUTHOR(S): Ismeshko, A. V.; Bowdrill, T.; Belosludtsev, Y. Y.; Hagan, M.  
 CORPORATE SOURCE: Baylor College of Medicine, Houston, TX, 77030, USA  
 SOURCE: Nucleic Acids Research (2001), 29(14), 3051-3058

Searched by Barb O'Bryen, STIC 308-4291

PUBLISHER: CODEN: NARHAD; ISSN: 0047-1347  
DOCUMENT TYPE: Oxford University Press  
JOURNAL: Journal  
LANGUAGE: English

AB The double helix is known to form as a result of hybridization of complementary nucleic acid strands in aq. soln. In the helix the neg. charged phosphate groups of each nucleic acid strand are distributed helically on the outside of the duplex and are available for interaction with cationic groups. Cation-coated glass surfaces are now widely used in biotecnol., esp. for covalent attachment of cDNAs and oligonucleotides as surface-bound probes or microarrays. These cationic surfaces can bind the nucleic acid backbone electrostatically through the phosphate moiety. Here we describe a simple method to fabricate DNA microarrays based upon adsorptive rather than covalent attachment of oligonucleotides to a pos. charged surface. We show that such adsorbed oligonucleotide probes form a densely packed monolayer, which retains capacity for base pair-specific hybridization with a soln. state DNA target strand to form the duplex. However, both strand disson. kinetics and the rate of DNase digestion suggest, on symmetry grounds, that the target DNA binds to such adsorbed oligonucleotides to form a highly asym. and unwound duplex. Thus, it is suggested that, at least on a charged surface, a non-helical DNA duplex can be the preferred structural isomer under std. biochem. conditions.

BT 13822-56-5, 3-Aminopropyltrimethoxysilane  
RL: AKG (Analytical reagent use); ANST (Analytical study); USES  
(Uses

(oligonucleotides form duplex with non-helical properties on pos.  
charged surface)

EN 13822-56-5 CAPLUS

UN 1-Propanamine, 3-(trimethoxysilyl)- (9CI) (CA INDEX NAME)

OMe

OMe Si (CH<sub>3</sub>)<sub>3</sub> NH<sub>2</sub>

OMe

REFERENCE COUNT: 11 THERE ARE 11 CITED REFERENCES AVAILABLE FOR THIS  
RECORD. ALL CITATIONS AVAILABLE IN THE EE FORMAT

14. ANSWER 13 OF 41 CAPLUS COPYRIGHT 2002 ACS

ACCESSION NUMBER: 2001:447194 CAPLUS

DOCUMENT NUMBER: 138173337

TITLE: A novel biosensor of DNA immobilization on nano-gold  
modified ITO for the determination of mifepristone  
AUTHOR(S): Xu, Jinchong; Zhu, Jun-Jie; Zhu, Yanling; Gu, Kai;  
Chen, Hong-Yuan

ABSTRACT SOURCE: Department of Chemistry, State Key Laboratory of  
Coordination Chemistry, Nanjing University, Nanjing,  
210003, Peop. Rep. China

SOURCE: Analytical Letters (2001), 34(4), 503-512

CODEN: ANALBP; ISSN: 0003-2719

PUBLISHER: Marcel Dekker, Inc.

DOCUMENT TYPE: Journal

LANGUAGE: English

AB A novel DNA modified indium tin oxide (ITO) electrode has been prepd.  
(3-Aminopropyl)Trimethoxysilane, gold nano-particles and DNA mols. are  
modified on the ITO electrode surface by self-assembly and electrochem.  
techniques, resp. This is a simple, stable, repeatable approach. The  
modified electrode can be used to detect mifepristone. A linear  
dependence of the peak currents on mifepristone concns. is obsd. in the  
range 4.times.10<sup>-7</sup>-6.times.10<sup>-6</sup>mol/L. The relative std. deviation is 4.5.



for six successive detns. at 1.times.10-6 mol/L soln. The detection limit is 1.times.10-7 mol/L.

IT 13822-56-5, (3-Aminopropyl)trimethoxysilane  
FL: ARU (Analytical role, unclassified); DEV (Device component use);  
ANST (Analytical study); USES (Uses).

(DNA immobilization on laser-gold modified ITO for detn. of mifepristone)

RN 13822-56-5 CAPLUS

CN 1-Propylamine, 3-(trimethoxysilyl) - (SCI) (CA INDEX NAME)

ONE

MeG Si (CH<sub>3</sub>)<sub>3</sub> NH<sub>2</sub>

ONE

REFERENCE COUNT: 10 THREE ARE 16 CITED REFERENCES AVAILABLE FOR THIS RECORD. ALL CITATIONS AVAILABLE IN THE RE FORMAT

I45 ANSWER 14 OF 41 CACHED OCTOBER 1992 ACS

ACCESSION NUMREF: 2001:10037 CAPLUS

DOCUMENT NUMBER: 14:10037

TITLE: A factorial analysis of silanization conditions for the immobilization of oligonucleotides on glass surfaces

AUTHOR(S): Halliwell, Catherine M.; Cass, Anthony E. G.  
CORPORATE SOURCE: Department of Biochemistry Imperial College of Science Technology and Medicine, University of London, London, SW7 2AY, UK

SOURCE: Analytical Chemistry (2001), 73(11), 2476-2483  
CITER: ANCHAM; ISSR: 0007-1700

PUBLISHER: American Chemical Society

DOCUMENT TYPE: Journal

LANGUAGE: English

AB The modification of glass surfaces with (3-aminopropyl)trimethoxysilane and the application of this to DNA chip technol. are described. A range of factors influencing the silanization method, and hence the no. of surface-bound, chem. active thiol groups, were investigated using a design of expt. approach based on anal. of variance. The no. of thiol groups introduced on glass substrates were measured directly using a specific radiolabel, [14C] cysteamine hydrochloride. For liq.-phase silanization, the no. of surface-bound thiol groups was found to be dependent on both postsilanization thermal curing and silanization time and relatively independent of silane concn., reaction temp., and sample pretreatment. Depending on the conditions used in liq.-phase silanization, (1.3-9.0) .times. 10<sup>12</sup> thiol groups/cm<sup>2</sup> on the glass samples were bound. The reliability and repeatability of liq.- and vacuum-phase silanization were also investigated. Eighteen-base oligonucleotide probes were covalently attached to the modified surfaces via a 3'-amino modification on the DNA and subsequent reaction with the crosslinking reagent N-(gamma-maleimidebutyryl)xy succinimide ester (GMBSE). The resulting probe levels were detd. and found to be stoichiometric with that of the introduced thiol groups. These results demonstrate that silanization of glass surfaces under specific conditions, prior to probe attachment, is of great importance in the development of DNA chips that use the simple concept of the covalent attachment of pre-synthesized oligonucleotides to silicon oxide surfaces.

IT 919-30-2, (3-Aminopropyl)trimethoxysilane

FL: ARU (Analytical role, unclassified); DEV (Device component use);  
ANST (Analytical study); USES (Uses)

(factorial anal. of silanization conditions for immobilization of

1. Nucleotides on glass surfaces!  
IN 919-30-2 CAPLUS  
UN 1-Propylamine, 3-(triethoxysilyl)- (9CI) (CA INDEX NAME)

CE+

EN 91-174,175 NH2

CE+

REFERENCE COUNT: 42 THERE ARE 42 CITED REFERENCES AVAILABLE FOR THIS RECORD. ALL CITATIONS AVAILABLE IN THE PE FORMAT

LA ANSWER 15 OF 41 CAPLUS COPYRIGHT 2002 ACS

ACCESSION NUMBER: 2001:194413 CAPLUS

EXHIBIT NUMBER: 199:4:176

TITLE: Peptide and small molecule **microarray** for high throughput cell adhesion and functional assays

AUTHOR(S): Halsey, James R.; Senil, M.; Park, Steven; Li, Snijun; Lam, Kit S.

LABORATORY SOURCE: UC Davis Cancer Center Division of Hematology/Oncology and Department of Internal Medicine, University of California Davis, Sacramento, CA, 95817, USA

SOURCE: Biocconjugate Chemistry (2001), 12(3), 346-353

CODEN: BOCHES; ISSN: 1043-1802

PUBLISHER: American Chemical Society

DOCUMENT TYPE: Journal

LANGUAGE: English

AB A novel class of chem. microchips consisting of glass microscope slides was prepd. for the covalent attachment of small mol. ligands and peptides through site-specific oxime bond or thiazolidine ring ligation reaction. Com. available microscope slides were thoroughly cleaned and derivatized with (3-aminopropyl)triethoxysilane (APTES). The amino slides were then converted to glyoxylyl deriva. via two different routes: (1) coupling of Fmoc-Ger followed by deprotection and oxim., or (2) coupling with protected glyoxylic acid and final deprotection with HCl. Biotin or peptide ligands derivatized at the carboxyl terminus with a 4,7,10-trioxa-1,13-tridecanediamine succinimic acid linker and an amino-oxy group or a 1,2-amino-triol group (e.g., cysteine with a free N.alpha.-amino group) were printed onto these slides using a DNA microarray spotter. After chem. ligation, the microarray of immobilized ligands was analyzed with three different biol. assays: (1) protein-binding assay with fluorescence detection, (2) functional phosphorylation assay using [ $\gamma$ - $^{32}$ P]-ATP and specific protein kinase to label peptide substrate spots, and (3) adhesion assay with intact cells. In the cell adhesion assay, not only can we det. the binding specificity of the peptide against different cell lines, we can also det. functional cell signaling of attached cells using immunofluorescence techniques in situ on the microchip. This chem. microchip system enables us to rapidly analyze the functional properties of numerous ligands that we have identified from the "one-bead-one-compd." combinatorial library method.

IN 919-30-2, (3-Aminopropyl)triethoxysilane

RL: ARJ (Analytical role, unclassified); DEW (Device component use);

ANST (Analytical study); USES (Uses)

peptide and small mol. **microarray** for high throughput cell adhesion and functional assays

IN 919-30-2 CAPLUS

UN 1-Propylamine, 3-(triethoxysilyl)- (9CI) (CA INDEX NAME)

OEt

EtO Si (CH<sub>2</sub>)<sub>3</sub> NH<sub>2</sub>

OEt

REFERENCE COUNT: 42 THERE ARE 42 CITED REFERENCES AVAILABLE FOR THIS RECORD. ALL CITATIONS AVAILABLE IN THE RE FORMAT

L45 ANSWER 16 OF 41 CAPLUS COPYRIGHT 2002 ACS

ACCESSION NUMBER: 2001:159116 CAPLUS

DOCUMENT NUMBER: 154:307437

TITLE: Controlled immobilization of DNA molecules using chemical modification of mica surfaces for atomic force microscopy: Characterization in air

AUTHOR(S): Omemura, Kizuo; Ishikawa, Mitsuru; Kurada, Reiko

CORPORATE SOURCE: Joint Research Center for Atom Technology (JRCAT)-Angstrom Technology Partnership (ATP), Tsukuba, Ibaraki, 305-0846, Japan

SOURCE: Analytical Biochemistry (2001), 295(2), 232-237

CODEN: ANBCA2; ISSN: 0003-2687

PUBLISHER: Academic Press

DOCUMENT TYPE: Journal

LANGUAGE: English

AB Immobilization of biomols. on surfaces while keeping the max. conformational flexibility of the mols. is one of the most important techniques for at. force microscopy imaging. We have developed two methods of controlling adsorption of DNA mols. on mica surfaces. The first method is the use of a mica surface modified with dild. 3-aminopropyltriethoxysilane (APS). Here we named this a "dild. AP3-treated mica (AP-mica)" technique. The second method is the use of a mica surface modified with mixed self-assembled monolayers of organosilanes. In both of the techniques, the no. of DNA mols. immobilized on a mica surface was controlled. Further, a conformational change of circular DNA, from a supercoiled to a relaxed form was obsd. for the mols. immobilized on a dild. AP-mica surface, when 254-nm UV light was irradiated. This observation demonstrated that flexibility of circular DNA mols. was kept on a dild. AP-mica surface. (c) 2001 Academic Press.

17 919-30-2, 3-Aminopropyltriethoxysilane

PL: APL (Analytical role, unclassified); DEV (Device component use);

ANST (Analytical study); USES (Uses)

(DNA immobilization using chem. modification of mica surfaces for at. force microscopy: characterization in air)

EN 919-30-2 CAPLUS

CN 1-Propylamine, 3-(triethoxysilyl)- (SCI) (CA INDEX NAME)

OEt

EtO Si (CH<sub>2</sub>)<sub>3</sub> NH<sub>2</sub>

OEt

REFERENCE COUNT: 21 THERE ARE 21 CITED REFERENCES AVAILABLE FOR THIS RECORD. ALL CITATIONS AVAILABLE IN THE RE FORMAT

L45 ANSWER 17 OF 41 CAPLUS COPYRIGHT 2002 ACS

ACCESSION NUMBER: 2001:152167 CAPLUS

DOCUMENT NUMBER: 154:159906

TITLE: Method for the covalent immobilization and labeling of biopolymers especially the preparation of nucleic acids

microarrays  
 INVENTOR(S): Anselme, Wilhelm; Faulstich, Konrad  
 PATENT ASSIGNEE(S): Europäisches Laboratorium Fuer Molekularbiologie  
 (EMBL), Germany  
 PUBLICATION: ECT Int. Appl., 33 pp.  
 CODEN: PIMPD2  
 DOCUMENT TYPE: Patent  
 LANGUAGE: German  
 FAMILY APP. NUM. COUNT: 1  
 PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
WO 2000-1014585	A1	20000301	WO 2000-EP8193	20000822

W: AE, AG, AL, AM, AT, AU, AV, BA, BF, BG, BH, BY, BZ, CA, CH, CN,  
 CR, CU, CZ, DE, DK, DU, EE, EG, FI, GE, GD, GE, GH, GM, HR,  
 HU, IB, IL, IN, IS, IT, KE, KG, KH, KR, LC, LK, LR, LS, LT,  
 LU, LV, MA, MD, MG, MF, MH, MW, MX, MY, NZ, OC, NZ, PL, PT, RO, RU,  
 SE, SE, SG, SI, SK, SL, TC, TM, TR, TT, TS, UA, UG, US, VE, VN,  
 YU, ZA, ZW, AM, AZ, BY, EG, EZ, ME, RU, TC, TM

EW: GR, SN, HE, LS, MW, ME, SD, SL, SZ, TG, UG, ZW, AT, BE, CH, CY,  
 DE, DK, ES, FI, FR, GE, GR, IE, IT, LU, MC, NL, PT, SE, BF, BJ,  
 CF, CG, CI, CM, GA, GN, GW, ML, ME, NE, SN, TD, TG

DE 12016073	A1	20010301	DE 2000-10016073	20000831
-------------	----	----------	------------------	----------

EP 1212466	A1	20020612	EP 2000-200356	20000822
------------	----	----------	----------------	----------

R: AT, BE, CH, DE, DK, EG, FR, GB, GE, IT, LI, LU, NL, SE, MC, PT,  
 IE, SI, LT, LV, FI, RO, MK, CY, AL

PRIORITY APPL. INFO.:

DE 1999-19946377 A 19990824

DE 1001-10016073 A 20000831

WO 2000-EP8193 W 20000822

AB The invention relates to methods for covalent immobilization of  
 biopolymers, esp. those of nucleic acids, on a solid phase. Covalent  
 bonds are made between primary or and secondary amino groups of said  
 biopolymers and groups of the solid phase which react with said amino  
 groups. Silica-based solid phases with defined functional groups are used  
 for the immobilization of 5' amino-modified nucleotides; the prepd. DNA  
 microarrays are used in amplification procedures.

IT 51895-58-0

RL: DEV (Device component use); USES (Uses)

(method for covalent immobilization and labeling of biopolymers esp.  
 prepn. of nucleic acid **microarrays**)

EN 51895-58-0 CAPLUS

CI 1,6-Hexanediamine, N-[3-(trimethoxysilyl)propyl]- (9CI) (CA INDEX NAME)

OMG

MOO S1 (CH2)3 NH (CH2)6 NH2

OMG

REFERENCE COUNT:

6

THERE ARE 6 CITED REFERENCES AVAILABLE FOR THIS  
 RECORD. ALL CITATIONS AVAILABLE IN THE IE FORMAT

LI ANSWER 19 OF 41 CAPLUS COPYRIGHT 2001 AFS

PUBLICATION NUMBER: 2001:3114 0 CAPLUS

PUBLICATION NUMBER: 194:70366

TITLE: Oligonucleotide arrays for high resolution HLA typing

INVENTOR(S): Petersdorf, Efi; W.; Guo, Zhen; Hansen, John A.;

Hood, Leroy

PATENT ASSIGNEE(S): Fred Hutchinson Cancer Research Center, USA;

University of Washington

Searched by Barb O'Bryen, STIC 308-4291

SOURCE: PCT Int. Appl., 83 pp.  
CODEN: PIXXD2  
DOCUMENT TYPE: Patent  
LANGUAGE: English  
FAMILY ACC. NUM. COUNT: 1  
PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
WO 2000/019036	A1	20001228	WO 2000-US16722	20000616

W: AU, CA, CH, US

RW: AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL,  
PT, SE

PRIORITY APPLIC. INFO.: US 1999-139843P P 19990617

AB Arrays of HLA Class I oligonucleotide probes on a solid support are provided, wherein the probes are sufficient to represent at least 80% of the known polymorphisms in exons 2 and 3 of the HLA Class I locus.

II 13822-56-5, Aminopropyltrimethoxysilane

EL: APF (Analytical role, unclassified); BUJ (Biological use, unclassified); THU (Therapeutic use); ANST (Analytical study); BIOL (Biological study); USES (Uses)

(solid support derivatized with; oligonucleotide

arrays for high resolu. HLA typing and transplant compatibility anal.)

EN 13822-56-5 CAPLUS

CN 1-Propanamine, 3-(trimethoxysilyl)- (PCI) (CA INDEX NAME)

OMe

MeO Si (CH<sub>2</sub>)<sub>3</sub> NH<sub>2</sub>

OMe

REFERENCE COUNT: 3 THERE ARE 3 CITED REFERENCES AVAILABLE FOR THIS RECORD. ALL CITATIONS AVAILABLE IN THE RE FORMAT

I46 ANSWER 1 OF 41 CAPLUS COPYRIGHT 2002 ACS

ACCESSION NUMBER: 2000:6989:9 CAPLUS

DOCUMENT NUMBER: 134:232351

TITLE: Covalent attachment of DNA to glass supports using a new silane coupling agent and chemiluminescent detection

AUTHOR(S): Zhang, Guojun; Zhou, Yikai; Wu, Xiaoyan; Yuan, Jinwei; Ren, Shu

CORPORATE SOURCE: Institute of Environmental Medicine, Tongji Medical University, Wuhan, 430030, Peop. Rep. China

SOURCE: Journal of Tongji Medical University (L030), 20(2), 89-91

CODEN: JTMUEI; ISSN: 0257-716X

PUBLISHER: Tongji Medical University

DOCUMENT TYPE: Journal

LANGUAGE: English

AB A new kind of silane coupling agent, N-(.beta.-aminoethyl)-.gamma.-aminopropyl triethoxysilane, was used for DNA direct attachment on the surfaces of glass supports, then the immobilized DNA was hybridized with horseradish peroxidase (HRP)-labeled probe, and detected by using enhanced chemiluminescent method. In comparison with .gamma.-aminopropyl triethoxysilane, the detection limits (S/N) of DNA were 10 pg and 75 pg resp. Several exptl. conditions of DNA attaching to glass supports were investigated, and the system of hybridization of nucleic acid on the surfaces of glass supports was developed.

IT 919-30-2, 3-APTES

RL: ANST (Analytical role, unclassified); BAC (Biological activity or effect), except adverse; BPR (Biological process); BSU (Biological study, unclassified); ANST (Analytical study); BIOL (Biological study); PROC (Process)

covalent attachment of DNA to glass supports using a new silane

coupling agent and chemiluminescent detection

RI 319-8-2 CAPLUS

TI 1-Propanamine, 3-(triethoxysilyl)- (9CI) (CA INDEX NAME)

OEt

EM: Si (CH<sub>3</sub>)<sub>3</sub> NH<sub>2</sub>

Et

REFERENCE COUNT: 5 THERE ARE 5 CITED REFERENCES AVAILABLE FOR THIS RECORD. ALL CITATIONS AVAILABLE IN THE RE FORMAT

44. ANSWER 20 OF 41 CAPLUS COPYRIGHT 2002 ACS

ACCESSION NUMBER: 2000:679668 CAPLUS

DOCUMENT NUMBER: 134:159699

TITLE: Protein microarrays for monitoring of structural changes of proteins via surface enhanced metal nano cluster resonance

AUTHOR(S): Mayer, Christian; Palkovits, Roland; Bauer, Georg; Schalkhammer, Thomas

CORPORATE SOURCE: Kluwer C. for Biotechnology, TU-Delft, Delft, 2628BC, Neth.

SOURCE: Micro Total Analysis Systems 2000, Proceedings of the 4th. MTAS Symposium, 4th, Enschede, Netherlands, May 14-18, 2000 (2000), 551-556. Editor(s): Van den Berg, Albert; Olthuis, W.; Bergveld, Piet. Kluwer Academic Publishers: Dordrecht, Neth.

CODEN: 69AJFE

DOCUMENT TYPE: Conference

LANGUAGE: English

AB Structural changes of ultra thin protein layers caused by changes in microenvironment, meaning a conformational change of the protein, were transduced into a optical signal obsd. directly as a color change of a biochip surface. We have successfully coated proteins as thin films of 10 to 100 nm onto optically reflecting ultra-flat and ultra-pure chip-surfaces via microdotting, spin-coating and subsequent photocrosslinking. The optical resonance effect was obtained by deposition of metal nano-clusters on top of the proteins. The response of this protein biochip array was measured spectroscopically in the visible and IR range of the spectrum. This set-up enabled us to transduce a change of protein conformation of various serum proteins and enzymes into a signal quant., reversibly and directly visible to the human eye.

11 3179-76-8

RL: NWU (Other use, unclassified); USES (Uses)

protein **microarrays** for monitoring of structural changes of proteins via surface enhanced metal nano cluster resonance)

RI 3179-76-8 CAPLUS

TI 1-Propanamine, 3-(diethoxymethylsilyl)- (9CI) (CA INDEX NAME)

(E)

PATENT INFORMATION:

CN 2(1H)-Pyrimidinone, 4-amino- (9CI) (CA INDEX NAME)

||  
|| NH<sub>2</sub>

BN 75-41-1 CAPLUS

CI 6H-Purin-6-one, 2-amino-1,7-dihydro- (9CI) (CA INDEX NAME)

BN ||  
N N

N NH

273752-55-9DP, immobilized 273752-56-0DP,  
immobilized 273752-57-1DP, immobilized  
273752-58-2DP, immobilized 273752-59-3DP,  
immobilized 273752-60-6DP, immobilized  
273752-61-7DP, immobilized 273752-62-8DP,  
immobilized 273752-63-9DP, immobilized


EL: DEV (Device component use); PEP (Physical, engineering or chemical  
process); ECT (Reactant); SPN (Synthetic preparation); PREP (Preparation);  
PROC (Process); RACT (Reactant or reagent); USES (Uses)  
(prepr. and detachment of; methods and comps. for performing arrays of  
chem. reactions on support surfaces using photoresists)

BN 273752-55-9 CAPLUS

CI L-Proline, N-[(1S)-1-carboxy-2-phenylethyl]-L-alanyl-,  
2-[1,1-dimethyl-3-[4-[2-oxo-2-[[3-(triethoxysilyl)propyl]amino]ethoxy]phen-  
yl]propyl] ester (9CI) (CA INDEX NAME)

Absolute stereochemistry.

PAGE 1-A






PAGE 1-B

EtO  
OEt  
Si  
OEt

RN 273752-56-0 CAPLUS  
CN L-Proline, N-[(1S)-1-carboxy-2-(2-nitrophenyl)ethyl]-L-alanyl-,  
2-[1,1-dimethyl-3-[4-[2-oxo-2-[[3-(triethoxysilyl)propyl]amino]ethoxy]phen-  
yl]propyl] ester (9CI) (CA INDEX NAME)

Absolute stereochemistry.

PAGE 1-A



PAGE 1-B


EtO  
OEt  
Si  
(CH<sub>2</sub>)<sub>3</sub> OEt

RN 273752-57-1 CAPLUS  
CN L-Proline, N-[(1S)-1,3-dicarboxypropyl]-L-alanyl-, 2-[1,1-dimethyl-3-[4-[2-  
oxo-2-[[3-(triethoxysilyl)propyl]amino]ethoxy]phenyl]propyl] ester (9CI)


(CA INDEX NAME)

Absolute stereochemistry.

PAGE 1-A




PAGE 1-B




273754-58-2 CAPLUS  
 L-Proline, N2-[(1S)-1-carboxy-2-phenylethyl]-N-(triphenylmethyl)-L-  
 asparaginyl-, 2-[1,1-dimethyl-3-[4-[2-oxo-2-[[3-  
 (triethoxysilyl)propyl]amino]ethoxy]phenyl]propyl] ester (9CI) (CA INDEX  
 NAME)

Absolute stereochemistry.

PAGE 1-A




PAGE 1-B



RN 273752-53-3 CAPLUS  
 CN L-Proline, N2-[(1S)-1-carboxy-2-(2-nitrophenyl)ethyl]-N-(triphenylmethyl)-  
 L-asparaginyl-, 2-[1,1-dimethyl-3-[4-[2-oxo-2-[[3-(triethoxysilyl)propyl]amino]ethoxy]phenyl]propyl] ester (9CI) (CA INDEX  
 NAME)

Absolute stereochemistry.

PAGE 1-A




PAGE 1-B

EtO  
OEt  
Si  
(CH<sub>2</sub>)<sub>3</sub> OEt


IN 273751-60-6 CAPLUS  
 IN L-Proline, N2-[(1S)-1,3-dicarboxypropyl]-N-(triphenylmethyl)-L-asparaginyl-  
 , 2-[1,1-dimethyl-3-[4-[2-oxo-2-[[3-(triethoxysilyl)propyl]amino]ethoxy]ph  
 enyl]propyl] ester (9CI) (CA INDEX NAME)

Absolute stereochemistry.

PAGE 1-A



PAGE 1-B



RN 273752-61-7 CAPLUS  
 CN L-Proline, N-[(1S)-1-carboxy-2-phenylethyl]-O-(1,1-dimethylethyl)-L-seryl-  
 , 1-[1,1-dimethyl-3-[4-[2-oxo-2-[[3-(triethoxysilyl)propyl]amino]ethoxy]ph  
 enyl]propyl] ester (9CI) (CA INDEX NAME)

Absolute stereochemistry.

PAGE 1-A




PAGE 1-B

EtO  
OEt  
Si  
OEt

RM 273751-62-8 CAPLUS  
 TI L-Proline, N-[(1S)-1-carboxy-2-(2-nitrophenyl)ethyl]-O-(1,1-dimethylethyl)-  
 L-seryl-, 2-[1,1-dimethyl-3-[4-[2-oxo-2-[[3-(triethoxysilyl)propyl]amino]e  
 thoxy[phenyl]propyl] ester (9CI) (CA INDEX NAME)

Absolute stereochemistry.

PAGE 1-A




PAGE 1-B

EtO  
OEt  
Si  
(CH<sub>2</sub>)<sub>3</sub> OEt

RN 273752-63-9 CAPLUS  
CN L-Proline, N-[(1S)-1,3-dicarboxypropyl]-O-(1,1-dimethylethyl)-L-seryl-,  
2-[1,1-dimethyl-3-[4-[2-oxo-2-[[3-(triethoxysilyl)propyl]amino]ethoxy]phen-  
yl]propyl] ester (9CI) (CA INDEX NAME)

Absolute stereochemistry.

PAGE 1-A



PAGE 1-B

EtO  
OEt  
CH<sub>2</sub>)<sup>2</sup> OEt

14 ANSWER 22 OF 41 CAPLUS COPYRIGHT 2002 ACS  
ACCESSION NUMBER: 2000:54088 CAPLUS  
DOCUMENT NUMBER: 132:90351  
TITLE: Photoluminescent semiconductor materials  
INVENTOR(S): Armstrong, David W.; Lafrance, Martine L.  
PATENT ASSIGNEE(S): Iatroquest Corporation, Can.  
SOURCE: ECT Int. Appl., 37 pp.  
CODEN: PIXXDE  
DOCUMENT TYPE: Patent  
LANGUAGE: English  
FAMILY ACC. NUM. COUNT: 1  
PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
WO 200006251	A1	20000110	WO 1999-CA642	19990709
W: AE, AL, AN, AT, AU, BA, BE, BG, BR, BY, CA, CH, CN, CU, CZ, DE, DK, EN, ES, FI, GB, GR, HE, HU, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LF, LS, LT, LU, LV, MD, MG, MK, MN, MW, MX, NO, NZ, PL, PT, RO, RU, SD, SE, SG, SI, SK, SL, TJ, TH, TR, TT, UA, UG, US, VN, YU, ZA, ZW, AM, AZ, BY, KG, KZ, MD, EU, TC, TM				
RW: GH, GM, KE, LS, MW, SD, SL, SS, UG, ZW, AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LI, MC, NL, PT, SE, BF, BJ, CF, CG, CI, CM, GA, GN, GW, ML, ME, NE, SN, TD, TG				
AU 9947682	A1	20000201	AU 1999-47682	19990709
EP 1198703	A1	20020414	EP 1999-210919	19990709
R: AT, BE, CH, DE, DK, EN, FR, GB, GR, IT, LI, LU, NL, SE, MC, PT, IE, FI, CY				

PRIORITY APPLN. INFO.:  
US 1991-00414P P 19980710  
WO 1999-CA642 W 19990709

AB Semiconductor materials having a porous texture are described which are modified with a recognition element and produce a photoluminescent response on exposure to electromagnetic radiation. The semiconductor materials may be doped, and they may be supported on a core material. The recognition elements, which can be selected from biomol., org., and inorg. moieties, interact with a target analyte to produce a modulated photoluminescent response, as compared with that of semiconductor materials modified with a recognition element only. The target analyte may be an inorg. or org. compd. or biomol., or an organism or a material derived from or produced by an organism. Methods for detecting an analyte are also described which entail comparing photoluminescence from the materials in a sample to that from the materials in the absence of a sample.



IT 919-30-2DP, .gamma.-Aminopropyltriethoxysilane, reaction products with oxidized porous silicon and recognition moieties 2530-83-8DP, 3-Glycidioxypropyltrimethoxysilane, reaction products with oxidized porous silicon and recognition moieties  
 EL: AEG (Analytical reagent use); SPN (Synthetic preparation); ANST (Analytical study); PREP (Preparation); USES (Uses)  
 opacoluminescent indicators based on surface-modified porous semiconductors)

RI 919-30-2 CAPLUS

CI 1-Propanamine, 3-(triethoxysilyl)- (9CI) (CA INDEX NAME)

OEt

E.O. Si (CH<sub>2</sub>)<sub>3</sub> NH<sub>2</sub>

OEt

RI 919-30-8 CAPLUS

CI Silane, trimethoxy[3-(oxiranylmethoxy)propyl]- (9CI) (CA INDEX NAME)

O

OMe

CH<sub>2</sub> C (CH<sub>2</sub>)<sub>3</sub> Si OMe

OMe

REFERENCE COUNT: 6 THERE ARE 6 CITED REFERENCES AVAILABLE FOR THIS RECORD. ALL CITATIONS AVAILABLE IN THE RE FORMAT

L15 ANSWER 01 OF 41 CAPLUS COPYRIGHT 2002 ACS

ACCESSION NUMBER: 1999:718211 CAPLUS

DOCUMENT NUMBER: 111:33271

TITLE: Chemically modified nucleic acids having enhanced lability towards solid supports, and uses thereof in high-density microarrays

INVENTOR(S): Bradley, Allan; Cai, Wei Wen

PATENT ASSIGNEE(S): Baylor College of Medicine, USA

SOURCE: PCT Int. Appl., 38 pp.

CODEN: PINXD2

DOCUMENT TYPE: Patent

LANGUAGE: English

FAMILY ACC. NUM. COUNT: 3

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
WO 9957323	A1	19991111	WO 1999-US9610	19990504
W: AU, CA, JP				
RW: AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE				
US 6048695	A	20000411	US 1998-71876	19980504
CA 2326684	AA	19991111	CA 1999-2326654	19990504
AU 9957861	A1	19991113	AU 1999-37361	19990504
EP 175544	A1	20010214	EP 1999-920342	19990504
R: AT, BE, CH, DE, DK, ES, FR, GB, GR, IT, LI, LU, NL, SE, MC, PT, IE, FI				
JP 2012513814	T2	20020514	JP 2000-547274	19990504
PRIORITY APPLN. INFO.:				
		US 1998-71876	A	19980504
		WO 1999-US9610	W	19990504

OTHER SOURCE(S): MARPAT 131:332971

AB The invention relates to novel chem. modified nucleic acids with enhanced lability towards solid supports, such as glass. These modified nucleic acids can be readily affixed to solid supports, for instance, a glass surface, without first derivatising the glass surface. In certain embodiments, the chem. modified nucleic acids of the invention are so modified via (1) compds. having a ring ether and an alkoxysilane group, (2) compds. having an amino group and an alkoxysilane group, (3) halogenated silanes, or (4) amine-contg. silanes reacted with brominated nucleic acids. High-d. microarrays based on these modified nucleic acids as well as methods for prepg. these microarrays are also useful.

BT 919-30-2DP, 3-Aminopropyltriethoxysilane, bound to a nucleic acid  
2530-83-8DP, 3-Glycidoxypropyltrimethoxysilane, bound to a nucleic acid

EL: ARG (Analytical reagent use); BPN (Biosynthetic preparation);  
ANST (Analytical study); BIOL (Biological study); PREP (Preparation); USES (Uses)

(chem. modified nucleic acids having enhanced lability towards solid supports, and uses thereof in high-d. microarrays)

EN 919-30-2 CAPLUS

TI 1-Propanamine, 3-(triethoxysilyl)- (9CI) (CA INDEX NAME)

OEt

BT SI (CH<sub>2</sub>)<sub>3</sub> NH<sub>2</sub>

OEt

EN 2530-83-8 CAPLUS

TI Silane, trimethoxy[3-(oxiranylethoxy)propyl]- (9CI) (CA INDEX NAME)

OMe

CH<sub>2</sub> O (CH<sub>2</sub>)<sub>3</sub> Si OMe

OMe

BT 71-30-7, Cytosine

EL: BOC (Biological occurrence); BSU (Biological study, unclassified);  
BIOL (Biological study); OCCU (Occurrence)

(modified nucleic acid comprising; chem. modified nucleic acids having enhanced lability towards solid supports, and uses thereof in high-d. microarrays)

EN 71-30-7 CAPLUS

TI 2(1H)-Pyrimidinone, 4-amino- (9CI) (CA INDEX NAME)

H  
N

BT 1591-21-5 14867-28-8, 3-Iodopropyltrimethoxysilane  
70892-80-7, 3-Bromooctyltrichlorosilane 82985-34-0,  
3-Bromooctyltrimethoxysilane

EL: ANU (Analytical role, unclassified); BUU (Biological use,

unclassified); ANST (Analytical study); BIOL (Biological study);  
USES (Uses)

Use in modifying nucleic acids; chem. modified nucleic acids having  
enhanced lability towards solid supports, and uses  
thereof in high-d. microarrays)

EN 1591-11-5 CAPLUS

CH Silane, dichloro(4-chlorobutyl)methyl- (7CI, 8CI, 9CI) (CA INDEX NAME)

Cl

Me Si (CH<sub>2</sub>)<sub>4</sub> Cl

Cl

EN 14867-28-8 CAPLUS

CH Silane, (3-iodopropyl)trimethoxy- (7CI, 8CI, 9CI) (CA INDEX NAME)

OMe

MeO Si (CH<sub>2</sub>)<sub>3</sub> I

OMe

EN 7-612-80-7 CAPLUS

CH Silane, (6-bromooctyl)trichloro- (9CI) (CA INDEX NAME)

Cl

Cl Si (CH<sub>2</sub>)<sub>3</sub> Br

Cl

EN 8-985-14-0 CAPLUS

CH Silane, (8-bromooctyl)trimethoxy- (9CI) (CA INDEX NAME)

OMe

MeO Si (CH<sub>2</sub>)<sub>8</sub> Br

OMe

REFERENCE COUNT: 2 THERE ARE 2 CITED REFERENCES AVAILABLE FOR THIS  
RECORD. ALL CITATIONS AVAILABLE IN THE RE FORMAT

145 ANSWER 24 OF 41 CAPLUS COPYRIGHT 2002 ACS

ACCESSION NUMBER: 1999:619530 CAPLUS

DOCUMENT NUMBER: 131:297347

TITLE: Addressable protein arrays on solid supports using  
capture oligonucleotides and RNA-protein fusions

INVENTOR(S): Kaimelis, Robert G.; Wagner, Richard

PATENT ASSIGNEE(S): Paylos, Inc., USA

SOURCE: PCT Int. Appl., 57 pp.

CODEN: PIXXD2

DOCUMENT TYPE: Patent

LANGUAGE: English

Searched by Barb O'Bryen, STIC 308-4291

FAMILY APP. NUM. COUNT: 1  
PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
WO 9911773	A1	19991014	WO 1999-037203	19990331
W: AL, AM, AT, AU, AS, BA, BB, BG, BF, BY, CA, CH, CN, CU, CZ, DE, DK, EE, ES, FI, GE, GL, GR, GU, HD, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MD, MG, MK, MN, MW, MX, NO, NZ, PL, PT, RG, RU, SD, SE, SG, SI, SK, SL, TC, TM, TR, TT, UA, UG, UK, VN, YU, ZA, ZW, AM, AZ, BY, EG, GT, HE, HU, IE, FI				
RW: GR, GM, HE, IG, MW, ND, NL, NZ, OG, SW, AT, BE, CH, CY, DE, DK, ES, FI, FR, GE, GR, HE, IL, LU, MC, ML, PL, SE, SF, SJ, SG, CC, CI, CM, GA, GN, GW, ML, MS, NE, SN, TD, TG				
CA 2213635	AA	19991014	CA 1999-2321635	19990331
AP 9934636	A1	19991025	AP 1999-34636	19990331
EP 1999-34636	A1	23-19117	EP 1999-34636	19990331
R: AT, BE, CH, DE, DK, ES, FR, GE, GR, IT, LI, LU, NL, SE, MC, PT, IE, FI				
JP 2002510505	T2	23020409	JP 2000-541454	19990331
PRIORITY APPLN. INFO.:				
US 1998-30686E P 19980403				
WO 1999-037203 W 19990331				
AB Disclosed herein are arrays of nucleic acid-protein fusions which are immobilized to a solid surface through capture probes which include a non-nucleosidic spacer group and an oligonucleotide sequence to which the fusion (such as an RNA-protein fusion) is bound. RNA-protein fusions are synthesized by in vitro translation of mRNA pools contg. a peptide acceptor such as puromycin attached to their 3'-ends, such that a covalent amid bond forms between the 3'-end of the mRNA and the C-terminus of the protein which it encodes. The arrays are prepd. by fixing oligonucleotide sequences, the capture probes, to a support in a defined array; the capture probes are then used to bind nucleic acid-protein fusions through base pairing between the nucleic acid component of the fusion and a complementary capture probe. The result of the binding interactions between the fusions and the capture probes is a defined, addressable array of proteins attached to a solid support. Also disclosed herein are solid supports on which these arrays are immobilized as well as methods for their prepn. and use (for example, for screening for protein-compd. interactions such as protein-therapeutic compd. interactions). Exemplary fusion chips are generated for FLAG, Ha11, and c-Myc epitope fusions.				
13822-56-5				
RI: DEV (Device component use); ECT (Reactant); EACT (Reactant or reagent); USES (Uses)				
(addressable protein arrays on solid supports using capture oligonucleotides and RNA-protein fusions)				
RN 13822-56-5 CAPLUS				
CI 1-Propanamine, 3-(trimethoxysilyl)- (9CI) (CA INDEX NAME)				

OMe

1999 SI (CH<sub>3</sub>)<sub>3</sub> NH<sub>2</sub>

OMe

REFERENCE COUNT: 7 THERE ARE 7 CITED REFERENCES AVAILABLE FOR THIS RECORD. ALL CITATIONS AVAILABLE IN THE RE FORMAT

14 ANSWER 25 OF 41 CAPLUS COPYRIGHT 2002 ACS  
SESSION NUMBER: 1999:458939 CAPLUS  
DOCUMENT NUMBER: 141:164790

Searched by Barb O'Bryen, STIC 308-4291


TITLE: Preparation and evaluation of p-tert-butylcalix[4]arene-bonded silica stationary phases for high-performance liquid chromatography  
AUTHOR(S): Xiao, Xiang-Zhu; Feng, Yu-Qi; Da, Shi-Lu; Zhang, Yan  
CORPORATE SOURCE: Dep. Chemistry, Wuhan Univ., Wuhan, 430072, Peop. Rep. China  
SOURCE: Chromatographia (1999), 49(11/12), 643-648  
CODEN: CHRGB7; ISSN: 0009-5893  
PUBLISHER: Friedrich Vieweg & Sohn Verlagsgesellschaft mbH  
DOCUMENT TYPE: Journal  
LANGUAGE: English

AB A method is proposed for prepn. of a 4-tert-butylcalix[4]arene-bonded silica stationary phase. Chem. modified 4-tert-butylcalix[4]arene is attached to silica gel by using [ $\gamma$ -(ethylenediamino)propyl]triethoxysilane as coupling reagent. The bonded phase was characterized by  $^{29}\text{Si}$  and  $^{13}\text{C}$  cross polarization/magic angle spinning solid-state NMR. The retention behavior of polycyclic arom. hydrocarbons (PAHs), nucleosides, and nucleobases was investigated on the bonded phase in the reversed-phase mode.

IT 71-30-7, Cytosine  
FL: ANT (Analyte); ANST (Analytical study)  
(prepn. and evaluation of tert-butylcalixarene-bonded silica stationary phases for HPLC)

EN 71-30-7 CAPLUS

CN 2-(1H)-Pyrimidinone, 4-amino- (9CI) (CA INDEX NAME)



N

IT 30858-91-4DP, [ $\gamma$ -(Ethylenediamino)propyl]triethoxysilane, reaction product with silica gel and tert-butyl[(chlorocarbonyl)methoxy]hydroxy-calixarene

FL: ARU (Analytical role, Unclassified); SPN (Synthetic preparation);

ANST (Analytical study); PREP (Preparation)

(prepn. and evaluation of tert-butylcalixarene-bonded silica stationary phases for HPLC)

EN 30858-91-4 CAPLUS

CN 1,2-Ethanediamine, N,N'-bis[3-(triethoxysilyl propyl)]- (9CI) (CA INDEX NAME)

OEt

OEt

EtO Si (CH<sub>2</sub>)<sub>3</sub> NH CH<sub>2</sub> CH<sub>2</sub> NH (CH<sub>2</sub>)<sub>3</sub> Si OEt

OEt

OEt

REFERENCE COUNT: 29 THERE ARE 29 CITED REFERENCES AVAILABLE FOR THIS RECORD. ALL CITATIONS AVAILABLE IN THE RE FORMAT

L43 ANSWER 16 OF 41 CAPLUS COPYRIGHT 2002 ACE

ACCESSION NUMBER: 1999:428814 CAPLUS

DOCUMENT NUMBER: 131:211144

TITLE: Atomic force microscopy imaging of DNA covalently immobilized on a functionalized mica substrate

AUTHOR(S): Shlyakhtenko, Luda S.; Gall, Alexander A.; Weimer, Jeffrey J.; Hawn, David D.; Lyubchenko, Yuri L.

Searched by Barb O'Bryen, STIC 308-4291

DEPARTMENT SOURCE: Department of Microbiology, Arizona State University,  
Tempe, AZ, 85287-2701, USA  
SOURCE: Biophysical Journal (1999), 77(1), 568-576  
CODEN: BIOJAU; ISSN: 0006-3495  
PUBLISHER: Biophysical Society  
DOCUMENT TYPE: Journal  
LANGUAGE: English

AB A procedure for covalent binding of DNA to a functionalized mica substrate is described. The approach is based on photochem. crosslinking of DNA to immobilized psoralen derivs. A tetrafluorophenyl (TFP) ester of tri-Na psoralen (trioxalen) was synthesized, and the procedure to immobilize it onto a functionalized aminopropyl mica surface (AP-mica) was developed. DNA mols. were cross-linked to trioxalen moieties by UV irradiation of complexes. The steps of the sample prepn. procedure were analyzed with XPS (XIS). Results from XPS show that an AP-mica surface can be formed by vapor phase deposition of silane and that this surface can be derivatized with trioxalen. The derivatized surface is capable of binding of DNA mols. such that, after UV crosslinking, they withstand a thorough rinsing with PBS. Observations with at. force microscopy showed that derivatized surfaces remain smooth, so DNA mols. are easily visualized. Linear and circular DNA mols. were photochem. immobilized on the surface. The mols. are distributed over the surface uniformly, indicating rather even modification of AP-mica with trioxalen. Generally, the shapes of sup-molled mols. electrostatically immobilized on AP-mica and those photocross-linked on trioxalen-functionalized surfaces remain quite similar. This suggests that UV crosslinking does not induce formation of a noticeable no. of single-stranded breaks in DNA mols.

BT 919-30-2

EL: AR (Analytical role, unclassified); ANST (Analytical study)  
(mica surface coated with,; imaging of DNA by at. force microscopy  
based on covalent photochem. crosslinking of DNA to trioxalen  
immobilized onto mica surface)

KN 919-30-3 CAPLUS

UN 1-Propanamine, 3-(triethoxysilyl)- (PCI) (CA INDEX NAME)

OE+

BT 919-30-3 NH2

OE+

REFERENCE COUNT: 44 THESE ARE 44 CITED REFERENCES AVAILABLE FOR THIS  
RECORD. ALL CITATIONS AVAILABLE IN THE RE FORMAT

INFO ANSWER 27 OF 41 CAPLUS COPYRIGHT 2002 ACS

ACCESSION NUMBER: 1999:11258 CAPLUS

DOCUMENT NUMBER: 130:143546

TITLE: Novel methods of attaching probes to a solid support  
and uses thereof

INVENTOR(S): Okamoto, Tadashi; Yamamoto, Nobuko; Suzuki, Tomochiro

PATENT APPLICANT(S): Canon Kabushiki Kaisha, Japan

SOURCE: Eur. Pat. Appl., 43 pp.

CODEN: EPKXNDW

DOCUMENT TYPE: Patent

LANGUAGE: English

FAMILY ACCL. NUM. COUNT: 1

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
EP 0 600 982	A2	19990208	EP 1998-336107	19980741

EP 345082 A3 19990311

E: AT, BE, CH, DE, DK, ES, FR, GE, GR, IT, LI, LU, NL, SE, MC, PT,  
IE, SI, LT, LV, FI, RO

JP 11117900 A1 19990718

JP 1998-209923 19980724

JP 2001066305 A1 20010316

JP 2000-232206 19980724

PRIORITY APPLN. INFO.:

JP 1997-207837 A 19970801

JP 1997-237046 A 19971020

JP 1998-209923 A 19980724

OTHER SOURCE(S): MARPAT 130:149846

AB Provided is a method of attaching probes to a solid support in a markedly high d. and efficiency. An extremely small amt. of probe is contained within a liq., and droplets of the liq. are delivered to the solid support via an ink jet ejection method, thereby forming spots which contain the probe. Since one or more substances can bind specifically to target probes and said probes are arranged in a large no. on a solid support, the method can be used to swiftly and accurately det. a base sequence of a nucleic acid or detect a target nucleic acid in a sample.

IT 1760-24-3, KBM603 2530-83-8, KBM403

EL: ECT (Reactant); EACT (Reactant or reagent)

Several methods of attaching probes to a **solid support**  
and uses thereof.

EN 1760-24-3 CAPLUS

CN 1,2-Ethanediamine, N-[3-(trimethoxysilyl)propyl]- (9CI) (CA INDEX NAME)

OMe

MeO Si (CH<sub>2</sub>)<sub>3</sub> NH CH<sub>2</sub> CH<sub>2</sub> NH<sub>2</sub>

OMe

EN 2130-11-8 CAPLUS

CN Silane, trimethoxy[3-(oxiranylmethoxy)propyl]- (9CI) (CA INDEX NAME)

O

OMe

CH<sub>2</sub> O (CH<sub>2</sub>)<sub>3</sub> Si OMe

OMe

145 ANSWER 28 OF 41 CAPLUS COPYRIGHT 2002 ACS

ACCESSION NUMBER: 1999:06139 CAPLUS

DOCUMENT NUMBER: 199:163949

TITLE: Novel polyethylenimine-based biomolecule arrays

INVENTOR(S): Van Ness, Jeffrey; Tabone, John C.; Moynihan, Kristen

PATENT ASSIGNEE(S): RapiGene, Inc., USA

SOURCE: PCT Int. Appl., 50 pp.

CODEN: PIXXD2

DOCUMENT TYPE: Patent

LANGUAGE: English

FAMILY ACC. NUM. COUNT: 1

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
WO 9904896	A1	19990204	WO 1998-US13246	19980721
W: AL, AM, AT, AU, BA, BB, BG, BR, BY, CA, CH, CN, CU, CZ, DE, DK, EE, ES, FI, GB, GE, GR, HU, IL, IS, JP, KE, KG, KI, KR, KZ, LC, LR, LE, LS, LT, LU, LV, MD, MG, MK, MN, NW, MX, NC, NZ, PL, PT,				

BO, BU, SD, SE, SG, SI, SK, TJ, TM, TR, TT, UA, UG, US, UZ, VN,  
 YU, AM, AZ, BY, BG, BZ, BU, BU, TM  
 BW: GH, GN, KE, LS, MW, SD, SZ, UG, ZW, AT, BE, CH, CY, DE, DK, ES,  
 FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE, BF, BJ, CF, CG, CI,  
 CM, GA, GN, GW, ML, MR, NE, SN, TD, TG

AT 88-5825 A1 19990116 AU 1998-95925 19980721

AT 88-5825 B2 20010712

EP 88-5825 A1 20000117 EP 1998-95925 19980721

B: AT, BE, CH, DE, DK, ES, FR, GB, GR, IE, IT, LU, NL, SE, MC, PT,  
 IE, FI

US 81-119 A 20001121 US 1998-120386 19980721

JP 80-151072 T2 20010707 JP 1998-50343 19980721

SECURITY AVELN. INFO.: US 1997-533528 P 19970722

WO 1998-051246 W 19980721

AB An array of biomols. is formed from a flat solid substrate, whereby said surface is covered with a layer of polyethylenimine (PEI) and this layer is divided among a plurality of discrete first regions abutted and surrounded by a contiguous second region. The process includes the step of depositing a biomol. into the first regions while maintaining the second region substantially free of the biomol.

BT 2530-83-8, 3-(2,3-Epoxypropoxy)propyltrimethoxysilane  
 EL: AN (Analytical role, unclassified); RCT (Reactant); ANST (Analytical study); RACT (Reactant or reagent)

use as bifunctional coupling agent; novel polyethylenimine-based biomol. arrays

BT 2530-83-8 CAPLUS

BT Silane, trimethoxy[3-(oxiranymethoxy)propyl]- (PCI) (CA INDEX NAME)

OMe

CH<sub>2</sub> O (CH<sub>2</sub>)<sub>3</sub> Si OMe

OMe

REFERENCE COUNT: 2 THERE ARE 2 CITED REFERENCES AVAILABLE FOR THIS RECORD. ALL CITATIONS AVAILABLE IN THE EE FORMAT

BT ANSWER 29 OF 41 CAPLUS COPYRIGHT 2002 ACT

ACCESSION NUMBER: 1997:258651 CAPLUS

DOCUMENT NUMBER: 127:1389

TITLE: Covalent attachment of hybridizable oligonucleotides to glass supports

AUTHOR(S): Joos, Bodo; Kuster, Herbert; Cone, Richard

CORPORATE SOURCE: Div. Infectious Diseases, Univ. Hospital, Zurich, CH-8091, Switz.

SOURCE: Analytical Biochemistry (1997), 247(1), 96-101  
 CODEN: ANBCA2; ISSN: 0003-2697

PUBLISHER: Academic

DOCUMENT TYPE: Journal

LANGUAGE: English

AB A simple, rapid, and efficient method for the covalent binding of oligonucleotides to solid glass supports was developed. Glass slides were derivatized with aminophenyl or aminopropyl silanes and 3'-succinylated target oligonucleotides were attached by carbodiimide-mediated coupling. Approx. 40-50% of the applied target oligonucleotides covalently bound to the derivatized glass. Hybridizations with radioactively labeled oligonucleotide probes showed that up to 90% of the attached oligonucleotides were available for hybridization. This system can conveniently be applied for studies on hybridization and detection of nucleic acids.

BT 919-30-2DP, 3-Aminopropyltrimethoxysilane, reaction products with



glass slides

EL: ARU (Analytical role, unclassified); DEV (Device component use); SPN (Synthetic preparation); ANST (Analytical study); PREP (Preparation); USES (Uses)

(covalent attachment of hybridizable oligonucleotides to glass supports)

RN 918-A-2 CAPLUS

CN 1-Propanamine, 3-(trimethylsilyl)- (901) (CA INDEX NAME)

OE:

EtO Si (CH<sub>2</sub>)<sub>3</sub> NH<sub>2</sub>

OE:

L45 ANSWER 30 OF 41 CAPLUS COPYRIGHT 2002 ACS

ACCESSION NUMBER: 1996-057311 CAPLUS

DOCUMENT NUMBER: 123:00153

TITLE: Carbazine dyes and derivatives for pH measurement

INVENTOR S.: Smith, Roger E.

PATENT ASSIGNER(S): Utah Medical Products, Inc., USA

SOURCE: U.S., 23 pp.

CODEN: USKXAM

DOCUMENT TYPE: Patent

LANGUAGE: English

FAMILY APP. NUM. COUNT: 1

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
US 5567614	A	19961021	US 1996-429622	19960427
CA 2219117	AA	19961021	CA 1996-2219117	19960426
WO 96-4284	A1	19961031	WO 1996-US5777	19960426
W: AU, AM, AT, AU, BG, BR, BY, CA, CH, CN, DE, DK, EE, ES, FI, GB, GE, HU, IS, JP, KE, KG, KP, KR, KZ, LK, LR, LS, LT, LU, LV, MD, MG, MK, MN, MW, MX, NO, NZ, PL, PT, RO, RU, SD, SE, SG, SI				
RW: BR, BS, BW, SD, SI, UG, AT, BE, CH, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE, BF, BG, CF, CG, CI, CM, GA, GN				
AU 6645757	A1	19961113	AU 1996-55757	19960426
AS 66797	B2	19961022		
GB 2314606	A1	19960107	GB 1996-22470	19960426
GB 2314606	B2	19960916		
DE 1981163	T	19960401	DE 1996-19611363	19960426
DE 1981163	C1	20020711		

PRIORITY APPLN. INFO.: US 1996-429622 A 19960427  
WO 1996-US5777 W 19960426

AB A compn. for detg. pH of a soln. comprises a fluorescent carbazine dye covalently bound to a solid support. A method of detg. pH of a soln. comprises placing the compn. in the soln., contacting the compn. with a selected wavelength of light to excite fluorescence by the carbazine dye, measuring intensities of the fluorescence at two selected wavelengths, calcn. a ratio of fluorescence intensities at the two selected wavelengths, and correlating the ratio with a predetd. relation of such ratios to pH. A fiber optic system for measuring pH of a soln. with the carbazine-dye-contd. compn. is also disclosed.

IT 2530-83-8

EL: RCT (Reactant); RACT (Reactant or reagent)

(prepn. of carbazine dyes and derivs. bonded to solid supports for pH measurement)

2530-83-8 CAPLUS  
Silane, trimethoxy[3-(oxiranylmethoxy)propyl]- (9CI) (CA INDEX NAME)

OMe

CH<sub>3</sub> O (CH<sub>2</sub>)<sub>3</sub> Si OMe

OMe

ANSWER 31 OF 41 CAPLUS COPYRIGHT 2002 ACS  
ACCESSION NUMBER: 1991:461921 CAPLUS  
DOCUMENT NUMBER: 117:1921  
TITLE: Oligonucleotide hybridizations on glass supports: a novel linker for oligonucleotide synthesis and hybridization properties of oligonucleotides synthesized in situ

AUTHOR(S): Masdos, Uwe; Southern, Edwin M.  
CORPORATE SOURCE: Dep. Biochem., Univ. Oxford, Oxford, OX1 3QU, UK  
SOURCE: Nucleic Acids Res. (1992), 20(7), 1679-84  
CODEN: NARHAD; ISSN: 0305-1048

DOCUMENT TYPE: Journal  
LANGUAGE: English

AB: A novel linker for the synthesis of oligonucleotides on a glass support is described. Oligonucleotides synthesized on the support remain tethered to the support after ammonia treatment and are shown to take part in sequence-specific hybridization reactions. These hybridizations were carried out with oligonucleotides synthesized on ballottini solid sphere glass beads and microscope slides. The linker has a hexaethylene glycol spacer, bound to the glass via a glycidoxypropyl silane, terminating in a primary hydroxyl group that serves as starting point for automated or manual oligonucleotide synthesis.

2530-83-8

EL: USES (Uses)

(glass support immobilization of, reaction with diols after, for synthesis of **solid support**-bound linker for oligonucleotide synthesis)

2530-83-8 CAPLUS

Silane, trimethoxy[3-(oxiranylmethoxy)propyl]- (9CI) (CA INDEX NAME)

OMe

CH<sub>3</sub> O (CH<sub>2</sub>)<sub>3</sub> Si OMe

OMe

ANSWER 32 OF 41 CAPLUS COPYRIGHT 2002 ACS  
ACCESSION NUMBER: 1991:097102 CAPLUS  
DOCUMENT NUMBER: 116:207102  
TITLE: Thymine bonded-stationary phase for high performance liquid chromatography  
AUTHOR(S): Zhu, Tao; Wang, Qinwei; Shen, Lianzhu; Lu, Chengxun; Fan, Yiliang  
CORPORATE SOURCE: Dep. Chem., Peking Univ., Beijing, 100871, Peop. Rep. China  
SOURCE: Chin. Chem. Lett. (1991), 2(7), 543-6  
CODEN: CCLEB7  
DOCUMENT TYPE: Journal

LANGUAGE: English

AB A new type of HPLC stationary phase contg. thymine deriv. was successfully prepd. It was found to give selective sepn. of nucleic acid bases and several purine derivs., such as caffeine and theophylline. The retention behavior and elution order of the solutes were interpreted in terms of mol. structure.

IT **919-30-2DP**, reaction products with silica gel and subsequently with thymineylpropionic acid-hydroxynorbornenedicarboximide reaction product

EL: SYN (Synthetic preparation ; **ANST (Analytical study)**; **PREP (Preparation)**)

(prepn. and use of, as stationary phase for sepn. of nucleic acid bases)

EN 919-30-2 CAPLUS

CN 1-Propanamine, 3-(triethoxysilyl)- (9CI) (CA INDEX NAME)

OEt

EtO Si (CH<sub>2</sub>)<sub>3</sub> NH<sub>2</sub>

OEt

IT **73-40-5, Guanine**

EL: **ANST (Analytical study)**

(sepn. of, from nucleic acid bases by HPLC on thymine bonded silica gel)

EN 73-40-5 CAPLUS

CN 6H-Purin-6-one, 2-amino-1,7-dihydro- (9CI) (CA INDEX NAME)

H<sub>2</sub>N      H  
          N      N

N                  NH

O

L45 ANSWER 33 OF 41 CAPLUS COPYRIGHT 2002 ACS

ACCESSION NUMBER: 1989:590993 CAPLUS

DOCUMENT NUMBER: 111:190933

TITLE: Silica gel or metal oxide chromatographic material and its use

INVENTOR(F): Hammer, Richard Frederick

PATENT ASSIGNEE(S): Chromatochem, Inc., USA

SOURCE: Eur. Pat. Appl., 3: pp.

CODEN: EPXXDW

DOCUMENT TYPE: Patent

LANGUAGE: English

FAMILY ACC. NUM. COUNT: 1

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
EP 295073	A2	19881214	EP 1988-305217	19880608
EP 295073	A3	19901214		
EP 295073	B1	19970210		
E: AT, CH, DE, FR, GB, LI, NL, SE				
CA 1520718	A1	19950727	CA 1988-568784	19880607

JP 61-12054	A2	19890316	JP 1989-141451	19890608
JP 61-12054	B4	19890726		
AT 15,176	E	19870313	AT 1989-305217	19890608
US 4,146,601	A	19890331	US 1991-062393	19910402
US 6,112,774	B1	20011028	US 1993-261450	19930303

## PRIORITY APPLN. INFO.:

US 1,847-53988	A	19170608
US 1,883-137765	A	19190427
US 1,990-485866	B1	19900223
US 1,991-662393	A3	19910402
US 1,993-70154	B1	19930601
US 1,995-337414	B1	19950301
US 1,996-714523	B1	19960916
US 1,997-948448	B1	19971014

AB Chromatog. materials (SEX, SEXYL, and SEXY' [C = substantially noncompressible solid support; B = binding group; X = substantially nonionic hydrophilic spacer; Y = coupling group; Y' = activated coupling group; L = affinity ligand] are provided. The solid support is silica gel or other metal oxide or ceramic. A process for chromatog. sepn. and detection of a target substance with the title material is also provided. The chromatog. material is substantially free of nonspecific adsorption and is stable at high pH. PEG 600-propylsilica (4% .mu.m) was prepd. and activated with carbonyldiimidazole. The activated silica gel was reacted 1st with hydrazine, then with periodate-oxidized ovalbumin, and packed into a HPLC column. Serum from a rabbit immunized against ovalbumin was loaded onto the column. Following removal of nonbound serum components by washing, IgG was eluted with 2% HOAc contg. 0.1M NaCl. Identify of the eluted, purified IgG was confirmed by SDS-PAGE and Western blot anal.

BT 13883-39-1D, reaction products with silica gel

RL: ANST (Analytical study)

(in prepn. of stationary phase for affinity chromatog., pH stability in relation to)

BN 13883-39-1 CAPLUS

CN Silane, (3-bromopropyl)trichloro- (6CI, 8CI, 9CI) (CA INDEX NAME)

Cl

Cl Si (CH<sub>2</sub>)<sub>3</sub> Br

Cl

L45 ANSWER 34 OF 41 CAPLUS COPYRIGHT 2002 ACS

ACCESSION NUMBER: 1989:208929 CAPLUS

DOCUMENT NUMBER: 110:208929

TITLE: Manufacture of silanized hydroxyethyl methacrylate-ethylene glycol dimethacrylate copolymers and their use as solid supports for affinity chromatographic methods for use in medicine and pharmaceutical industry

INVENTOR(S): Schuessler, Werner; Coupek, Jiri; Hiepe, Falk

PATENT ASSIGNEE(S): Akademie der Wissenschaften der DDR, Ger. Dem. Rep.

ADDRESS: Ger. (East), 4 pp.

CODEN: GEXXA3

DOCUMENT TYPE: Patent

LANGUAGE: German

FAMILY APP. NUM. COUNT: 1

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
-----	----	-----	-----	-----

FD 256720 A1 19860518 FD 1986-286583 19860129

OTHER SOURCE(S): MARPAT 110:208929

AB A process for the manuf. of chem. activated hydroxyethyl methacrylate-ethylene glycol dimethacrylate copolymer (I) in the form of shaped objects comprises the treatment of I with organosilanes  $(XSi)_m nSiR_4-n$  (X = amino, CO, CO<sub>2</sub>, isothiocyano, epoxy, diazo, NCO, NO, sulfinydril, halocarbonyl; R1 = alkyl, alkylphenyl, Ph; R = alkyl, phenoxy, halo, m = 0-20, n = 1-3) and optionally with hetero- or monofunctional reagents. Macroporous I (Separon Hema-1000; particle size 15-25  $\mu$ m; inner surface 70 m<sup>2</sup>/g; mol. wt. exclusion 2 times 10<sup>6</sup>) (5 g) was incubated with 10% aminopropyltriethoxysilane (EB 1114) in 1:1 EtOH-H<sub>2</sub>O at pH 2.5 for 6 h at 60.degree., washed with EtOH-H<sub>2</sub>O and 0.1M phosphate buffer at pH 6.3, and the resulting gel was incubated with 5% glutaraldehyde for 2 h at 39.degree. and subsequently washed with phosphate buffer. The activated gel was incubated with human IgG (18.6 mg IgG/ml 0.1M phosphate buffer) for 2 h at 37.degree. and overnight at 4.degree.; 36.7 mg IgG/g (>95%) were bound on activated I.

IT 919-30-2DP, reaction products with Separon HEMA and glutaraldehyde  
2602-34-8DP, reaction products with Separon HEMA and (aminopropyl)triethoxysilane and glutaraldehyde  
EL: PPEP (Preparation)

quant. of, as **solid support** for affinity chromatog.)

RN 119-30-2 CAPLUS

CN 1-Propanamine, 3-(triethoxysilyl)- (9CI) (CA INDEX NAME)

OEt

EtO Si (CH<sub>2</sub>)<sub>3</sub> NH<sub>2</sub>

OEt

RN 2602-34-8 CAPLUS

CN Silane, triethoxy[3-(oxiranylmethoxy)propyl]- (9CI) (CA INDEX NAME)

O

OEt

CH<sub>2</sub> C (CH<sub>2</sub>)<sub>3</sub> Si OEt

OEt

L45 ANSWER 15 OF 41 CAPLUS COPYRIGHT 2002 ACS

ACCESSION NUMBER: 1989:150718 CAPLUS

DOCUMENT NUMBER: 110:150718

TITLE: Modification of silufol plate silica gel by amino groups of aminopropyltriethoxysilane and their use for separation of nucleic acids components

AUTHOR(S): Karpova, S. F.; Pupkova, V. I.; Khripin, Yu. L.

CORPORATE SOURCE: Sci.-Res. Design-Technol. Inst. Biol. Active Subst., Berdsk, USSR

SOURCE: Zh. Anal. Khim. (1989), 44(1), 127-30

CODEN: ZAKHA8; ISSN: 0044-4502

DOCUMENT TYPE: Journal

LANGUAGE: Russian

AB A simple method for modifying Silufol plate silica gel by amino groups of aminopropyltriethoxysilane involves submerging of the corn plates in 1-2% soln. of Me ethoxysilane in EtOH for 50-60 min. The plates are dried for 23-30 min at room temp. and washed once with EtOH. The sepn. selectivity

of these plates (for sugars, guanosine, and its phosphates) is not inferior when compared with Merck com. plates NH2-F254. Ribonucleotides, deoxyribonucleotides and impurities of nucleoside N bases and their phosphates were sep'd. by a mobile phase contg. AcOH and EtOH.

IT 73-40-5, Guanine 73-40-5D, Guanine,


nucleotides

RL: ANST (Analytical study)

sepn. of, by TLC, aminopropyltrimethoxysilane-modified silica gel for)

RN 73-40-5 CAPLUS

W 6H-Purin-6-one, 2-amino-1,7-dihydro- (9CI) (CA INDEX NAME)




N NH

O

RN 73-40-5 CAPLUS

W 6H-Purin-6-one, 2-amino-1,7-dihydro- (9CI) (CA INDEX NAME)



N NH

O

IT 919-30-2, Aminopropyltriethoxysilane

RL: ANST (Analytical study)

silica gel-modified with, for nucleic acid component sepn., by TLC)

RN 919-30-2 CAPLUS

W 1-Propanamine, 3-(triethoxysilyl)- (9CI) (CA INDEX NAME)

OEt

EtO Si (CH2)3 NH2

OEt

147 ANSWER 36 OF 41 CAPLUS COPYRIGHT 2002 ACS

ACCESSION NUMBER: 1989:17873 CAPLUS

DOCUMENT NUMBER: 110:17873

TITLE: Synthesis and characterization of highly stable bonded phases for high-performance liquid chromatography column packings

AUTHOR(S): Kirkland, J. J.; Glajch, J. L.; Farlee, R. P.

CORPORATE SOURCE: Exp. Sta., E. I. du Pont de Nemours and Co., Wilmington, DE, 19898, USA

JOURN: Anal. Chem. (1989), 61(1), 2-11

CODEN: ANCHAM; ISSN: 0003-2700

DOCUMENT TYPE: Journal

LANGUAGE: English

AB: Two new classes of silane-modified silicas were synthesized and

characterized by chromatog. and spectroscopic techniques. These new bonded phases are significantly more stable toward hydrolysis than previous bonded-phase silicas; retention and column efficiency are comparable. The first type uses bifunctional (or "bidentate") silanes contg. one reactive atom on each of two silicon atoms that connect through a bridging group such as -O- or -(CH<sub>2</sub>)<sub>n</sub>-. The second type uses a monofunctional silane with at least two bulky groups (e.g., isopropyl) on the silane silicon atom. These bulky groups provide steric protection to the Si-O-Si bond formed between the silane and the surface of the silica. The new bonded-phase silicas exhibit highly reproducible gradient elution high-performance sepns. of peptides and proteins with low-pH mobile phases.

IT 116698-58-9DP, reaction products with silica gels  
 117559-36-1DP, reaction products with silica gels  
 RI: ANST (Analytical study); PREP (Preparation)  
 (prepn. and characterization and evaluation of, as stationary phases in: HPLC for anal. with low-pH mobile phases)  
 RN 116698-58-9 CAPLUS  
 CN Silane, ethoxybis(1-methylethyl)[3-(oxiranylmethoxy)propyl]- (9CI) (CA INDEX NAME)

O OEt  
 CH<sub>2</sub> (CH<sub>2</sub>)<sub>3</sub> Si Pr-i  
 i-Pr

RN 117559-36-1 CAPLUS  
 CN 1-Propanamine, 3-[ethoxybis(1-methylethyl)silyl]- (9CI) (CA INDEX NAME)

OEt  
 i-Pr Si (CH<sub>2</sub>)<sub>3</sub> NH<sub>2</sub>  
 i-Pr

I45 ANSWER 37 OF 41 CAPLUS COPYRIGHT 2002 ACS  
 ACCESSION NUMBER: 1988:4034:0 CAPLUS  
 DOCUMENT NUMBER: 109:3469  
 TITLE: Multicoated ferromagnetic chromium dioxide particles for use as **solid support** in heterogeneous immunoassays and bioaffinity separations  
 INVENTOR(S): Lau, Hon Beng Phillip; Yang, Esther Koo; Jacobson, Edward Wayne  
 PATENT ASSIGNEE(S): du Pont de Nemours, E. I., and Co., USA  
 SOURCE: Eur. Pat. Appl., 26 pp.  
 CODEN: EPXNDW  
 DOCUMENT TYPE: Patent  
 LANGUAGE: English  
 FAMILY ACCL. NUM. COUNT: 1  
 PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
EP 246770	A1	19871014	EP 1987-103692	19870314
EP 246770	B1	19920112		
R: AT, BE, CH, DE, ES, FR, GB, GR, IT, LI, LU, NL, SE				
CA 1269673	A1	19911001	CA 1987-531885	19870312

AT 1987-103692	E	19920215	AT 1987-103692	19870314
ES 1987-103692	A1	19930211	ES 1987-103692	19870314
JP 1987-14117	A2	19871335	JP 1987-14117	19870314
JP 1987-14117	B4	19921013		
DK 8701367	A	19870919	DK 1987-1367	19870317
MAGNETIC AFFIN. INFO.:			US 1986-341107	19860318
			EP 1987-103692	19870314

AB CrO<sub>2</sub> particles are modified to have desirable characteristics as solid support materials for immunoassays or for bioaffinity sepsns. The particles are surface reduced and coated with protective silica and silane layers. Such treatment prevents hydrolytic degran. of the particles, and provides a functionalized coat. CrO<sub>2</sub> particles were surface reduced in an aq. soln. of NaHSO<sub>3</sub>, then treated with NaAlO<sub>2</sub> and Na<sub>2</sub>SiO<sub>3</sub> soln. contg. Na acetate, pH 9. The particles were coated with 3-aminopropyltriethoxysilane. The chromate leaching test of these particles gave an absorbance of 0.02 at 372 nm. The particle settling time was 8 min. In an immunoassay for the detn. of TSH, a serum sample was mixed with an enzyme-labeled anti-TSH .beta.-subunit monoclonal antibody (MAb), then mixed with a slurry of particles carrying anti-TSH .alpha.-subunit MABs. The immune complexes formed were removed magnetically. The complexes were resuspended in a substrate soln. and incubated, the absorbance of the quenched soln. was read. Human serum contg. 0, 5, 25, and 50 .mu.U TSH/mL gave an absorbance of 0.1135, 0.1829, 0.4836, and 0.7947, resp.

199-30-2, 3-Aminopropyltriethoxysilane 5089-72-5

RL: ANST (Analytical study)

Surface-reduced magnetic chromium dioxide particles coated with silica and, for immunoassays and bioaffinity sepsns.)

919-30-2 CAPLUS

1-Propanamine, 3-(triethoxysilyl)- (9CI) (CA INDEX NAME)

DEt

EtO Si (CH<sub>2</sub>)<sub>3</sub> NH<sub>2</sub>

DEt

9089-72-5 CAPLUS

1,2-Ethanediamine, N-[3-(triethoxysilyl)propyl]- (9CI) (CA INDEX NAME)

DEt

EtO Si (CH<sub>2</sub>)<sub>3</sub> NH CH<sub>2</sub> CH<sub>2</sub> NH<sub>2</sub>

DEt

11 ANSWER 38 OF 41 CAPLUS COPYRIGHT 2002 ACS

ABSTRACT NUMBER: 1993:403447 CAPLUS

INVENT NUMBER: 109:3447

TITLE: Analytical method and kit for detecting and measuring specifically sequenced nucleic acid using flubrescent intercalation compounds and waveguides as solid support

INVENTOR(S): Sutherland, Ronald Macdonald; Bromley, Peter; Gentile, Bernard

PATENT ASSIGNEE(S): Battelle Memorial Institute, Switz.

SOURCE: Eur. Pat. Appl., 30 pp.

CODEN: EPXNDW

Searched by Barb O'Brien, STIC 308-4291



DOCUMENT TYPE: Patent  
 LANGUAGE: English  
 FAMILY ACC. NUM. COUNT: 1  
 PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
EP 84100	A1	19871111	EP 1987-810004	19870418
RI: EE, GB, DE, FR, SE, IT, NL, BE				
WO 8700956	A1	19871119	WO 1987-EP234	19870502
W: AU, BR, DK, FI, JP, NO, US				
AU 8701034	A1	19871201	AU 1987-75838	19870502
JP 9110221	T2	19890126	JP 1987-503871	19870502
FI 8700176	A	19871210	FI 1987-5770	19871230
NO 8801011	A	19880210	NO 1988-10	19880104
DK 8800036	A	19880217	DK 1988-6	19880104
PRIORITY APPL. INFO.:			EP 1986-810201	19860505
			WO 1987-EP234	19870502

AB A waveguide coated with single-stranded probe nucleic acids and carrying an internally reflected wave signal is contacted with an analyte soln. contg. denatured test DNA or RNA and fluorescent marker dye. Analyte nucleic acid with sequences homologous to that of the probe polynucleotide will hybridize therewith with concomitant binding of the fluorescent dye to the resultant duplex structures. Fluorescence resulting from the interaction of the wave signal at the waveguide/analyte interface with the signal generating centers created within the space probed by the evanescent component of the wave signal is detected and provides useful information on said sequences homologous to that of the probe nucleic acids. A plate waveguide with poly(dA) affixed (prepn. described for oligo-dA on aminopropyl glass plate) was affixed into a flow cell and a base-line signal was obtained with buffer in the cell. Both ethidium bromide and poly-dT were mixed and injected into the flow cell and the reaction was monitored. In a control, only ethidium bromide was added. The monitoring reaction was effectively immediate and only specific interaction between double-stranded DNA was detected.

IT 919-30-2, 3-Aminopropyltriethoxysilane

EL: ANST (Analytical study)

(grafting of, on waveguide, for nucleic acid attachment, nucleic acid detn. in relation to)

EN 919-30-3 CAPLUS

CN 1-Propylamine, 3-(triethoxysilyl)- (9CI) (CA INDEX NAME)

Et

EtO Si (CH<sub>2</sub>)<sub>3</sub> NH<sub>2</sub>

Et

145 ANSWER 39 OF 41 CAPLUS COPYRIGHT 2002 ACS

ACCESSION NUMBER: 1988:34236 CAPLUS

DOCUMENT NUMBER: 108:34236

TITLE: Polymer-modified silica-based supports for gel permeation chromatography of biopolymers

AUTHOR(S): Ivanov, A. E.; Zigis, L.; Turchinskii, M. F.; Kop'ev, V. P.; Reshetov, R. D.; Zubov, V. P.; Kastrikina, L. N.; Lonskaya, N. I.

CORPORATE SOURCE: Inst. Biorg. Khim. im. Shemyakina, Moscow, USSR

SOURCE: Mol. Genet., Mikrobiol. Virusol. (1987), (11), 39-46, 1 plate  
 CODEN: MGEVDU

DOCUMENT TYPE: Journal  
LANGUAGE: Russian  
AB Macroporous glass treated with .gamma.-aminopropyltriethoxysilane and then with 1:1 copolymer of N-vinylpyrrolidone and acryloyl chloride was prepd. and used for sepn. of influenza, Sendai, etc. viruses. The sorbent possesses low absorption activity but had higher stability and better hydrodynamic properties than commonly used sorbents (Sephacrose 4B, porous glass). The sorbent can be used repeatedly without regeneration (>20 times) and could be regenerated by washing with 1:1 iso-PrOH-H2O, when the chromatog. properties are totally restored. The inert sorbent was also used for the sepn. of Escherichia coli tRNA from 70 S ribosomes.  
BT 919-30-2, .gamma.-Aminopropyltriethoxysilane  
RL: ANST (Analytical study)  
    Glass treatment with, copolymer modification after, for gel chromatog. support prepr.)  
BT 919-30-2 CAPLUS  
BT 1-Propanamine, 3-(triethoxysilyl)- (9CI) (CA INDEX NAME)

GE1

BT 31-38(1)-NH2

GE1

BT ANSWER 40 OF 41 CAPLUS COPYRIGHT 2002 ACS  
ACCESSION NUMBER: 1988:31015 CAPLUS  
DOCUMENT NUMBER: 188:31015  
TITLE: Alkoxy silanes for the preparation of silica based stationary phases with bonded polar functional groups  
AUTHOR(S): Engelhardt, Heinz; Orth, Peter  
CORPORATE SOURCE: Angew. Phys. Chem., Univ. Saarlandes, Saarbruecker, Fed. Rep. Ger.  
SOURCE: J. Liq. Chromatogr. (1987), 10(8-9), 1999-2022  
CODEN: JLCHDE; ISSN: 0143-1919  
DOCUMENT TYPE: Journal  
LANGUAGE: English  
AB For prepn. of polar bonded phases with alkoxy silanes, an activator and a catalyst are required to achieve surface coverages comparable to those obtained with chlorosilanes. For activation a monolayer of H2O on the silica surface is sufficient. The most active catalyst in many cases has been p-toluenesulfonic acid, however, for silanes with basic groups Et3N gives better coverages. Silanes with polar groups tend to adsorb also with this group onto the surface thus preventing chem. binding via alkoxy groups. Long time experiences in the prepn. of amino phases, anion and cation exchangers and hydrophilic bonded phases for protein anal. are summarized.  
BT 35141-36-7D, reaction products with silica  
RL: ANST (Analytical role, unclassified); ANST (Analytical study)  
    as stationary phase, for anion-exchange liq. chromatog.)  
BT 35141-36-7 CAPLUS  
BT 1-Propanaminium, N,N,N-trimethyl-3-(trimethoxysilyl)-, chloride (9CI) (CA INDEX NAME)

OMe

MeO Si (CH<sub>2</sub>)<sub>3</sub> N<sup>+</sup>Me<sub>3</sub>

OMe

● Cl<sup>-</sup>

IT 919-30-2D, 3-Aminopropyltriethoxysilane, reaction products with silica

EL: ARU (Analytical role, unclassified); ANST (Analytical study)  
(as stationary phases, for liq. chromatog.)

EN 919-30-2 CAPLUS

CN 1-Propanamine, 3-(triethoxysilyl)- (9CI) (CA INDEX NAME)

OEt

EtO Si (CH<sub>2</sub>)<sub>3</sub> NH<sub>2</sub>

OEt

IT 71-30-7, Cytosine 73-40-5, Guanine

EL: ANT (Analyte); ANST (Analytical study)  
(sepr. of, from nucleobases, chem.-bonded silica stationary phases for cation-exchange liq. chromatog.)

EN 71-30-7 CAPLUS

CN 2(1H)-Pyrimidinone, 4-amino- (9CI) (CA INDEX NAME)

O  $\begin{array}{c} \text{H} \\ \text{N} \end{array}$  NH<sub>2</sub>

N

EN 73-40-5 CAPLUS

CN 6H-Purin-6-one, 2-amino-1,7-dihydro- (9CI) (CA INDEX NAME)

E2N  $\begin{array}{c} \text{H} \\ \text{N} \end{array}$  N

N NH

O

L45 ANSWER 41 OF 41 CAPLUS COPYRIGHT 2002 ACS

ACCESSION NUMBER: 1986:417540 CAPLUS

DOCUMENT NUMBER: 105:17540

TITLE: Manipulation of stationary-phase acid-base properties  
by a surface-levelling effect. Boronic  
acid-saccharide complexation

AUTHOR(S): Lochmuller, C. H.; Hill, Walter B.

Searched by Barb O'Bryen, STIC 308-4291

ORIGINATOR: P. M. Gross Chem. Lab., Duke Univ., Durham, NC, 27706, USA  
JOURNAL: ACS Symp. Ser. (1986), 297 (Chromatogr. Sep. Chem.), 210-25  
CODEN: ACSMC8; ISSN: 0097-6156  
DOCUMENT TYPE: Journal  
LANGUAGE: English

AB The use of boronic acid-substituted, amine-modified silica gel stationary phases for the HPLC sepn. of saccharides and nucleosides under neutral conditions was studied. Five stationary phases were prepd. using Partisil 10. The capacity factors for selected saccharides and nucleosides on columns packed with these stationary phases are given. The presence of residual amine groups in the surface bound, silica-based phenylboronic acid phases lowers the apparent pKa of the acid groups. This surface buffering effect permits boronate-saccharide complexation to occur at much lower pH values than is typically the case.

BT 102712-18-5D, reaction products with silica gel

RL: ANST (Analytical study)

as stationary phases for high-performance liq. chromatog. sepn. of nucleosides and saccharides)

EN 102712-18-5 CAPLUS

BT Boronic acid, [4-[[[3-(ethoxydimethylsilyl)propyl]amino]methyl]phenyl]- (9CI) (CA INDEX NAME)

OEt

CH<sub>2</sub> NH (CH<sub>2</sub>)<sub>3</sub> Si Me

Me

HO B

SE

BT 73-40-5

RL: ANT (Analyte); ANST (Analytical study)

high-performance liq. chromatog. st. on boronic acid-substituted amine-modified silica gel stationary phases)

EN 73-40-5 CAPLUS

BT 6H-Purin-6-one, 2-amino-1,7-dihydro- (9CI) (CA INDEX NAME)

H  
H N N

N NH

919-30-2 18306-79-1

RL: RMT (Reactant); ANST (Analytical study)

reaction of, with silica gel)

EN 919-30-2 CAPLUS

BT 1-Propanamine, 3-(triethoxysilyl)- (9CI) (CA INDEX NAME)

OEt

EtO Si (CH<sub>2</sub>)<sub>3</sub> NH<sub>2</sub>

OEt

RN 18306-79-1 CAPLUS

CN 1-Propanamine, 3-(ethoxydimethylsilyl)- (9CI) (CA INDEX NAME)

OEt

Me Si (CH<sub>2</sub>)<sub>3</sub> NH<sub>2</sub>

Me

FILE 'HOME' ENTERED AT 14:58:20 ON 13 AUG 2002

